



KaVo SONICflex[®] Tips

KaVo SONICflex tips



KaVo. Dental Excellence.

888-ASK-KAVO · www.kavousa.com



SONICflex[®] handpiece overview




**MASTER
SERIES**

**EXPERT
SERIES**

Fast and comfortable
change of tips - only
one turn





Variety developed to a fine point.

The wide range of options offered by the SONICflex, will soon make it one of your most-used handpieces.

Choose the right tip for your needs, from the wide range of options.

Numerous possibilities – main features and advantages for you:

- Over 50 different tips
- For every indication
- Ensuring successful, minimally-invasive treatment

Your daily advantage in performance.

Your KaVo SONICflex is the most powerful airscaler of its class. The handpiece with a wide range of uses. It supports your efforts with adjustable performance levels, for every indication – strong and precise.


Superior technology – main features and advantages for you:

- Three performance settings – optimum working conditions, for every indication
- Substantially quieter than comparable models offered by competitors (up to 15 dB)
- No additional equipment is necessary; just connect it to the KaVo MULTIflex coupling, without an annoying, over-floor cables or additional foot controls
- Optimum illumination: 25,000 LUX – for clear visualization
- Legendary KaVo quality – the SONICflex 2008 L represents the third scaler generation from KaVo, the culmination of 25 years of knowledge and experience




KaVo. Dental Excellence.

Hygiene



- 4 tips: Universal, Sickle, Perio, Perio extra long (subgingival treatment)
- Universal application
- Gentle removal of calculus
- Pleasant patient treatment



SONICflex scaler tip no. 5

Treatment: for removal of large areas of calculus
Form: Universal



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.5171 (5)
1005.8949 (5A)

SONICflex scaler tip no. 6

Therapy: for removal of large areas of supragingival calculus
Form: Sickle




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.5181 (6)
1005.8950 (6A)

SONICflex scaler tip no. 7

Therapy: for removal of large areas of supragingival calculus
Form: Perio



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.5191 (7)
1005.8951 (7A)

SONICflex scaler tip no. 8

Therapy: for removal of supragingival and subgingival calculus
Form: Perio, extra-long





Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

0571.5371 (8)
1006.1953 (8A)

SONICflex scaler tip no. 48

Therapy: removal of plaque, cleansing of fissures, cleansing of cavities and crown preparations
Form: various attachments, single-use brushes




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L
clean brush no. 1
clean brush no. 1A
clean brush no. 2
clean brush no. 2A
clean brush no. 3
clean brush no. 3A
clean brush no. 4
clean brush no. 4A
clean brush no. 5
clean brush no. 5A
clean brush no. 6
clean brush no. 6A
tip no. 48 Clean Inset (w/o Brush)
tip no. 48 A Clean Inset (w/o Brush)

0571.0001
1006.2026
0571.0402 (1)
1004.4125 (1A)
0571.0412 (2)
1004.4126 (2A)
0571.0422 (3)
1004.4127 (3A)
0571.0432 (4)
1004.4128 (4A)
0571.0472 (5)
1004.4129 (5A)
0571.0482 (6)
1004.4130 (6A)
0571.0401 (7)
1006.1982 (7A)






Periodontology



- Protection of the collagenous soft tissue and the root surface, from injuries
- Effective cleansing in pockets up to 9 mm deep
- Improved bacterial control
- Effective removal of subgingival plaque and calculus around implants



SONICflex paro tip no. 60

Therapy: gentle pocket cleansing during periodontal treatment; effective cleansing in pockets up to 9 mm deep; for the anterior region, smooth surfaces
 Form: long straight tip - rounded surface



Order number: compatible with SONICflex 2003 / L
 compatible with SONICflex 2008 / L



0571.7402 (60)
 1006.1934 (60A)

SONICflex paro tip no. 61

Therapy: gentle pocket cleansing during periodontal treatment; effective cleansing in pockets up to 9 mm deep; from first/third quadrant to central incisor accessed buccally, then cleansing of front teeth, premolars and molars of the second/fourth quadrant, accessed palatally / lingually
 Form: angled tip, left - rounded surface




Order number: compatible with SONICflex 2003 / L
 compatible with SONICflex 2008 / L



0571.7412 (61)
 1006.1935 (61A)

SONICflex paro tip no. 62

Therapy: gentle pocket cleansing during periodontal treatment; effective cleansing in pockets up to 9 mm deep; from second/fourth quadrant to central incisor accessed buccally, then cleansing of front teeth, premolars and molars of the first/third quadrant, accessed palatally / lingually
 Form: angled tip, right - rounded surface



Order number: compatible with SONICflex 2003 / L
 compatible with SONICflex 2008 / L




0571.7422 (62)
 1006.1936 (62A)

SONICflex implant tip

Therapy: effective removal of subgingival plaque and calculus around implants
 Form: conical tips, single-use plastic base (PEEK)




Order number: compatible with SONICflex 2003 / L
 compatible with SONICflex 2008 / L
 single-use plastic base (30 pcs.)



1003.8167 (48)
 1006.2027 (48A)
 1003.8168



Minimally-Invasive Surgery




- Minimally-invasive therapy, for incipient caries
- Safe and ultra-conservative, caries excavation
- Maximum preservation of healthy tooth tissue
- Preparation from all angles possible
- Minimally-invasive preparation, in fissure patterns
- Patient loyalty, thanks to "less fear of drilling", perfect for the treatment of children

SONICflex micro tip no. 30

Therapy: minimally-invasive treatment, for incipient caries
Form: small hemisphere, **mesial**, diamond coating D46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L




0571.6811 (30)
1006.1969 (30A)

SONICflex micro tip no. 31

Therapy: minimally-invasive treatment, for incipient caries
Form: small hemisphere, **distal**, diamond coating D46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L




0571.6801 (31)
1006.1971 (31A)

SONICflex micro tip no. 32

Therapy: minimally-invasive treatment, for incipient caries
Form: large hemisphere, **mesial**, diamond coating D46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.6831 (32)
1006.1973 (32A)




SONICflex micro tip no. 33

Therapy: minimally-invasive treatment, for incipient caries

Form: large hemisphere, **distal**, diamond coating D46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.6821 (33)
1006.1975 (33A)


SONICflex bevel tip no. 58

Therapy: bevelling of approximal cavity edges

Form: **mesial**, diamond coating D46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.7392 (58)
1006.1988 (58A)

SONICflex bevel tip no. 59

Therapy: bevelling of approximal cavity edges

Form: **distal**, diamond coating D46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 L




0571.7382 (59)
1006.1990 (59A)

SONICflex seal tip no. 45

Treatment: minimally-invasive preparation of fissures, minimal discomfort when opening of fissures
Form: conical, straight tip 0.5 mm, diamond coating D 46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L




1000.8323 (45)
1007.1506 (45A)

SONICflex cariex D tip no. 42

Therapy: ultra-conservative, sonoabrasive caries excavation, fine
Form: ball, 0.8 mm, diamond coating D 46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L




1000.7105 (42)
1006.1980 (42A)

SONICflex cariex D tip no. 43

Therapy: ultra-conservative, sonoabrasive caries excavation, fine
Form: ball, 1.2 mm, diamond coating D 64



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



1000.7167 (43)
1006.1981 (43A)



SONICflex cariex TC tip no. 71

Therapy: controlled, layer-by-layer removal of carious tissues

Form: ball, 1.0 mm, hard metal



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

1000.7362 (71)
1006.2002 (71A)

SONICflex cariex TC tip no. 72

Therapy: controlled, layer-by-layer removal of carious tissues

Form: ball, 1.4 mm, hard metal




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L


1000.7363 (72)
1006.2004 (72A)



Prosthetics





- Precise preparation of various bevel angles, around margins
- Precise preparation result
- Reproducible results
- Improved preparation margins to assist dental technicians
- No damage to adjacent teeth



SONICflex prep gold tip no. 49, mesial

Therapy: precise preparation of various bevel-angles around margins, after "rough preparation" with rotary instrument; this enables improved dental restoration

Form: **mesial**, tips with defined bevel-angles: lateral - 60° cervical - 75°, diamond coating D 46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

0571.7212 (49)
1006.1983 (49A)

SONICflex prep gold tip no. 50 distal

Therapy: precise preparation of various bevel-angles around margins, after "rough preparation" with rotary instruments; this enables improved dental restoration

Form: **distal**, tips with defined bevel-edges: lateral - 60°, cervical - 75°, diamond coating D 46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

0571.7222 (50)
1006.1984 (50A)

SONICflex prep ceram tip no. 51

Treatment: precise preparation of various bevel-angles around margins, with gradual angle transitions from cervical to lateral areas

Form: **mesial**, tips with defined bevel-angles: lateral - 60°, cervical - 75°, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

0571.7262 (51)
1006.1985 (51A)




KaVo. Dental Excellence.


SONICflex prep ceram tip no. 52

Treatment: precise preparation of various bevel-angles around margins, with gradual angle transitions from cervical to lateral areas

Form: **distal**, tips with defined bevel-edges, lateral - 60°, cervical - 75°, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.7272 (52)
1006.1986 (52A)

SONICflex prep CAD / CAM tip no. 34

Therapy: perfect adaptation of cavity preparations, to the specific requirements of CAD-CAM systems

Form: **mesial**, lateral angle 60°, cervical angle 90°, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L




1002.1984 (34)
1006.1977 (34A)

SONICflex prep CAD / CAM tip no. 35

Therapy: perfect adaptation of cavity preparations, to the specific requirements of CAD-CAM systems

Form: **distal**, lateral angle 60°, cervical angle 90°, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L,



1002.1986 (35)
1006.1979 (35A)

SONICflex cem tip no. 12

Therapy: gentle and secure placement of inlays

Form: rounded plastic cover




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L
CEM attachment




0571.5431 (12)
1006.2021 (12A)
0571.7142

Root tip resection



- Correct axial preparation, of defined and medically exact retrograde cavities
- Protection of periapical bone
- Simplified access to modified apex, even with complex anatomical conditions



SONICflex retro tip no. 56

Therapy: searching for/locating of root canal; significant work simplification during retrograde cavity preparation

Form: extremely narrow cylindrical tip, **left**, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.7322 (56)
1006.2040 (56A)


SONICflex retro tip no. 57

Therapy: searching for/locating of root canal; significant work simplification during retrograde cavity preparation

Form: narrow cylindrical tip, **right**, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



0571.7332 (57)
1006.2041 (57A)


SONICflex retro tip no. 16

Therapy: preparation of root canal, significantly simplified access, even with highly-convoluted roots

Form: **left**, small cylinder, diamond coated D 64



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L




0571.5541 (16)
1006.2036 (16A)



SONICflex retro tip no. 17

Therapy: preparation of root canal, significantly simplified access, even with highly-convoluted roots
 Form: **right**, small cylinder, diamond coating D 64




Order number: compatible with SONICflex 2003 / L
 compatible with SONICflex 2008 / L

0571.5581 (17)
1006.2035 (17A)

SONICflex retro tip no. 20 and no. 21

Therapy: preparation of the undercut; correct, internal circumferential root canal preparation, with minimal bone defects
 Form: pear-shaped, **small, left** (no. 20), **right** (no. 21), diamond coating D 64



Order number: compatible with SONICflex 2003 / L
 matching for SONICflex 2008 / L

0571.5521 (20)
1006.2037 (20A)



Order number: compatible with SONICflex 2003 / L
 matching for SONICflex 2008 / L

0571.5561 (21)
1006.2038 (21A)

SONICflex retro tip no. 55

Therapy: for roots, with apices highly-inclined towards buccal cavity
 Form: diamond coating D 46




Order number: compatible with SONICflex 2003 / L
 compatible with SONICflex 2008 / L

0571.7342 (55)
1006.2039 (55A)



Endodontics




- Safe preparation of the canal entrance, with direct view when using a microscope
- Simplified, subsequent canal preparation
- Aligned tip set
- Safe and gentle removal of hard materials
- Improved location of canal entrances



SONICflex endo tip no. 66

Therapy: shaping of access cavity, removal of protuberances
Form: button, **large**, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L


1000.5825 (66)
1006.1992 (66A)

SONICflex endo tip no. 67

Therapy: expansion of larger channels
Form: conical **125°**, diamond coating D 46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



1000.5822 (67)
1006.1994 (67A)

SONICflex endo tip no. 68

Therapy: expansion of larger channels
Form: conical **112°**, diamond coating D 46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

1000.5823 (68)
1006.1996 (68A)



SONICflex endo tip no. 69

Therapy: removal of residual root canal filling materials
Form: button, **small**, diamond coating D 46




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

1000.5827 (69)
1006.1998 (69A)

SONICflex endo tip no. 70

Therapy: for particularly long and wide canals
Form: conical **117°**, diamond coating D 25 anterior area is not diamond coated




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L


1000.5821 (70)
1006.2000 (70A)



Surgical





- Minimizing of nerve and membrane injuries
- Gentle treatment, due to oscillating motion
- Highest safety with external sinus lift



SONICflex bone tip no. 80

Therapy: sectioning of Sharpey's fibers; extension of alveolae; removal of root residues; gentle extraction before implantation; bone splitting

Form: square, cutting



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L




1004.3875 (80)
1006.2006 (80A)


SONICflex bone tip no. 81

Therapy: preparation of bone window for sinus lift; with very massive bones, if applicable, rotating and sectioning of substantia corticalis, plus removal of remaining bone using the diamond-coated tip

Form: ball, large, diamond coating D 46



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



1004.3876 (81)
1006.2008 (81A)


SONICflex bone tip no. 83

Therapy: removal of bone blocks, osteotomies (distraction, bone splitting); if applicable, dealing with extremely thick, substantia corticalis

Form: saw sagittal



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L



1004.3878 (83)
1006.2012 (83A)



SONICflex bone tip no. 85

Therapy: the disc-shaped form, makes it ideal for safe repositioning of the Schneider membrane, during a sinus lift

Form: elephant foot




Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

1006.0645 (85)
1007.1624 (85A)

SONICflex bone tip no. 86

Therapy: ideal for easy and quick harvesting, of greater quantities of bone, especially useful for bone augmentation

Form: bone harvester





Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L

1006.0741 (86)
1007.1625 (86A)

SONICflex bone tip no. 87

Therapy: ideal for precise, extremely narrow incisions, when cutting ligaments

Form: saw axial, replaceable saw blades. Saw blade thickness: 0.15 mm



Order number: compatible with SONICflex 2003 / L
compatible with SONICflex 2008 / L
Saw blade (5 pcs.)

1006.1406 (87)
1007.1627 (87A)
1006.1405 (Set of 5)

Accessories:



Order number: **Torque Wrench - Tip Removal Tool**

1000.4887



Order number: **Sterilization Cassette**

0411.9101



Order number: **Key for Adapter**

0411.0892

Perio Attachment



Order number: coupling for the supply of therapeutic fluids such as sterile saline solution for SONICflex 2003 / L
coupling for the supply of therapeutic fluids such as sterile saline solution for SONICflex 2008 / L

1000.6296

1006.5966





KaVo Hygiene 360°

Everything Easy. Everything KaVo.

KaVo offers you the complete solution for every hygiene step:



KaVo SONICflex 2008 L



KaVo PROPHYflex 3



KaVo PROPHYwiz



DIAGNOdent pen



Hygienist Stool

Go to WWW.TRYKAVO.COM to start your free and easy trial today!

KaVo Dental
11727 Fruehauf Drive
Charlotte, NC 28273
888.ASK.KAVO (888.275.5286)
Fax 1.704.927.0816
www.kavousa.com

©2011 KaVo Dental Corporation. 906.9203/rev0/8.11



KaVo. Dental Excellence.