



Now you can hold, transport, store, refrigerate, incubate, and dispose of test tubes, all in one convenient unit.

TEST N' TOSS™ DISPOSABLE TEST TUBE RACK – NOW MADE FROM RECYCLED MATERIALS



The Test n' Toss™ disposable test tube rack is perfect for clinical labs, hospitals, research labs, reference labs, EMS, even physicians' offices and alternate care sites. The Test n' Toss, designed with the lab professional in mind, holds eighty 16 x 100 mm tubes in a leak-resistant plastic shell to contain spills and is inexpensive enough to dispose of when contaminated.

As an option, use a second shell as a cover, completely enclosing contents for transport or storage and allowing racks to be stacked. Uniquely designed divider shelves can be used singly or in pairs for extra tube stability.

- Size: Each rack holds eighty 16 x 100 mm test tubes
- Quantity: Carton of 10 half-shells and dividers, makes 10 open racks or 5 covered, stackable racks
- Dimensions: Open rack: 11"L x 9"W x 3"H
Covered rack: 11"L x 9"W x 5"H
- Shipping weight: 3 lbs.
- Reorder number: 18710



Whitney Medical Solutions
6153 C Mulford Street
Niles, IL 60714 • USA
(800) DEVICES • (800) 338-4237

Manage Your Lab
With WHITNEY MEDICAL SOLUTIONS:
Real Solutions for Real Problems™



Your Lab!
~~lab odors,~~
~~broken glass,~~
~~biohazards,~~

You'll wonder how you ever managed ¹ without us._{sm}

TERMINAL® AUTOCLAVE DEODORANT

Autoclave Deodorants

Shipping Weight: 2 lbs.
Quantity: Carton of 200 sachets
Order numbers:

Regular	BH3001
Fresh Lemon	BH3002
Mountain Pine	BH3003
Floral Bouquet	BH3004



Autoclave Deodorants eliminate lab odors and maintain a fresh lab environment with **FOUR FRESH SCENTS.**

SAFE-KEEPER™

You can save valuable space on your bench top with the hangable and leak-resistant Safe-Keeper™.



Safe-Keeper™

Dimensions: 5"L x 2"W x 11"H
Capacity: 1100cc
Shipping Weight: 13 lbs.
Quantity: Case of 200, With 1 Hook-Up™
Order number: BH2005
Features: Solid bottom construction avoids leaks; can be wall-mounted.

Hook-Up™

Dimensions: 2"L x 2"W x ¾"H
Shipping Weight: 1 lb.
Quantity: Box of 4
Order number: BH6012
Features: Self-adhesive plastic hook. Great for keys, light items, hanging extra Safe-Keepers™.

TERMINAL® STARTER SETS AND BIOHAZARD BAGS

2 Liter



7 Liter



18 Liter



60 Liter



Starter Sets

Implement your infectious waste management program with Terminal Biohazard and System Starter Sets. Sets consist of a reusable container & close-fitting lid with international biohazard symbol, 100 autoclavable bags sized to fit container, 100 closures, and instructions. Sets are available with red bags printed with the international biohazard symbol or unprinted blue bags.

Container	Bags	Unprinted Blue	Printed Red
2 liter	9"x17"	BH4201	BH4301
7 liter	14"x19"	BH4202	BH4302
18 liter	19"x23"	BH4203	BH4303
60 liter	25"x35"	BH4204	BH4304

Biohazard Bags

Eliminate risk of infection; protect staff & patients. Sterilize your infectious waste at temperatures up to 275°F (135°C). Treatment of infectious waste prior to disposal is required by many states and eliminates need for special disposal in some states. Available in red or blue with elastomeric closures. 9"x17" bags will fit most countertop sterilizers. Bags are sized to fit starter sets. 200 per case.

Size	Fits	Unprinted Blue	Printed Red
9"x17"	2 liter	BH1201	BH1301
14"x19"	7 liter	BH1202	BH1302
19"x23"	18 liter	BH1203	BH1303
25"x35"	60 liter	BH1204	BH1304

TERMINAL® INFECTIOUS WASTE CONTROL LOG



Documents the information required by OSHA and infectious waste disposal laws. Sample Infection Control Plan contains all you need to establish your biohazard management program: suggested procedures, training program, and employee certificates. Control Log, with 2-piece indicator tags, maintains a permanent record of treatment & disposal.

Log with 100 tags	BH7001
100 replacement tags	BH7002

For thirty years, we have designed, refined and manufactured convenient and cost-effective safety-based Whitney products for your lab environment.

Convenience - Where you need it, when you need it, with disposal management systems made for your workstation and common disposal areas.

Cost effective - Affordable, durable, disposable products designed to be tossed when you are finished.

Backed by on-site manufacturing, warehousing, distribution and responsive customer service, Whitney Medical Solutions meets the lab professional's needs from the smallest specialized lab to the larger lab environments of hospitals and industrial research.

As specialists in this area, we put our expertise to work for you.

Steve Whitney

Steven Whitney
President

BROKEN GLASS KEEPER™

Safely collect and dispose of non-biohazardous broken glass and sharp objects.



Floor Model

Dimensions: 12"W x 12"D x 27"H
Capacity: 63 liters
Shipping Weight: 14 lbs.
Quantity: Case of 6
Order number: BH2008
Features: Features wide-mouth opening to discard entire racks of glass tubes at once.

Bench-Top Model

Dimensions: 6"W x 6"L x 10"H
Capacity: 3 liters
Shipping Weight: 15 lbs.
Quantity: Case of 50
Order number: BH2004
Features: Lined with plastic to contain moderate amounts of liquid.



BIOHAZARD BENCH-TOP KEEPER™

Collect workstation waste economically without making unnecessary trips to your central biohazard.



Large (shown)

Dimensions: 6"L x 6"W x 10"H
Capacity: 3000cc
Shipping Weight: 15 lbs.
Quantity: Case of 50
Order number: BH2003
Features: Free standing; webbed corners and plastic lining keep liquids from soaking through.

Small

Dimensions: 4"L x 4"W x 7"H
Capacity: 1400cc
Shipping Weight: 12.5 lbs.
Quantity: Case of 100
Order number: BH2002
Features: Same great features as large size; takes up less bench space.

BIOHAZARD KEEPER™—FLOOR MODEL

All Terminal® brand Keepers™ are constructed of sturdy plastic-lined or corrugated paperboard, are latex free, autoclavable, and can be incinerated.

Dispose of non-sharp biohazardous lab waste including large items—even entire racks at once.



Terminal® brand Biohazard Keeper™

Dimensions: 12"W x 12"D x 27"H
Capacity: 63 liters
Shipping Weight: 14 lbs.
Quantity: Case of 6
Order number: BH2009
Features: Sturdy corrugated construction, heavy-duty 2 mil liner, and unique wide-mouth opening with Reclose-A/Remove-A-Lid™ cover.

BIOHAZARD PIPET KEEPER™

You can dispose of pipets and other long objects up to 15" in length easily with the Pipet Keeper™.



Pipet Keeper™

Dimensions: 3"W x 3"H x 15½"L
Capacity: 2250cc
Shipping Weight: 8 lbs.
Quantity: Case of 50
Order number: BH2006
Features: Collects contaminated long items right at your workstation.