



Digital Heat Gun

The variable temperature Digital Heat Gun offers 12 heat settings and automatic shut-off feature.



Features

- 12 temperature settings allows the user to choose the correct degree of heat needed for each splint making application.
- Cool down setting with automatic shut-off is ideal for safety at the splint table.
- Dial-in temperature with the digital touch-pad control panel.
- Heat gun adjusts from 250°F to 1350°F.
- Ergonomic handle, 1500 watts, 5100 BTU.
- Product Number P7530. *Product is yellow, not black.



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com