



3pp[®] Final Flexion Wrap[™]



The adjustable 3pp[®] Final Flexion Wrap[™] makes it easy to restore end range IP joint flexion. Soft foam-lining and wrap around design puts the user in control to make regaining full finger flexion comfortable and quick.


Features

- Hook receptive material allows the patient to adjust force to tolerance for fast, successful treatment.
- Soft foam lining distributes pressure and holds the wrap in place for improved control.
- As part of a hand rehabilitation program, this wrap is ideal for regaining full finger flexion.
- Washable. Latex-free.

Order Information

One size fits most.

Item Number	Size
P1004-5	Package of 5
P1004-25	Package of 25



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com