



3pp[®] Finger Trapper[™]



3pp[®] Finger Trappers[™] provides a secure hold for application of traction for post-op or dynamic splinting.

Features

- For use on fingers or thumb as part of a hand rehabilitation program.
- Puts an end to using fingernail hooks, tape or glue.
- Eyelet holds traction lines securely and makes it possible to apply distraction at the MP or PIP.
- Machine washable. Latex-free.

Order Information

One size fits most.

Item Number	Size
P1006-5	Package of 5
P1006-25	Package of 25

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com