



3pp[®] PF Lift[®]



The 3pp[®] PF Lift[®] relieves the pain of Plantar Fasciitis or heel spurs, day or night, with gentle controlled stretch. Straps apply adjustable force to lift and supinate the foot to relieve stress on the plantar fascia.

Features

- The soft, flexible 3pp PF Lift is designed to apply adjustable force across the bottom of the foot to hold a gentle stretch on the plantar fascia.
- Eases early morning pain.
- Can be worn in all styles of shoes and even when running.
- Non-slip foam-lined straps stay in place without sliding on the foot or leg.
- Hand or machine wash. Latex-free.
- For those with diabetes or poor circulation, it is recommended that wearers discuss this splint with a health care provider before wearing.

Order Information

Measure around the leg above the ankle. Fits right or left foot.

	Item Number
Sm/Med	P5006-23
Lrg/X-Lrg	P5006-45

Watch our **instructional video on the PF Lift** on our website:
www.3pointproducts.com/product-videos-of-splints/3pp-pf-lift-video



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com