



# ThumSaver™ CMC



The ThumSaver™ CMC Short immobilizes the thumb and stabilizes the wrist to reduce the pain of arthritis, tendinitis or overuse syndromes. This splint can also be used for deQuervain's tenosynovitis and Gamekeeper's thumb.

## Features

- Molded with thin, lightweight plastic, the ThumSaver CMC Short is a seamless splint that supports the thumb and limits wrist motion to reduce pain experienced when gripping or pinching.
- Carefully designed to support the thumb and wrist to allow protected hand function.
- Recommended for those experiencing thumb tendinitis and basal joint arthritis.
- Felt pads are included to provide cushioning where needed.

## Order Information

Measure the length from the IP joint (knuckle at the end of the thumb) to the wrist. Measure around the IP joint.

	<b>Right</b>	<b>Left</b>	<b>Length</b>	<b>Circumference</b>
Sm	P2000-R2	P2000-L2	3" to 3 1/2"	2 1/2" to 2 3/4"
Med	P2000-R3	P2000-L3	3 1/2" to 4"	2 3/4" to 3 1/8"
Lrg	P2000-R4	P2000-L4	4" to 4 1/2"	3 1/8" to 3 1/4"
X-Lrg	P2000-R5	P2000-L5	4 1/2" to 5"	3 1/4" to 3 1/2"



**Call 888-378-7763**

Mon-Fri 8:30am - 5:00pm EST for ordering details

**[www.3pointproducts.com](http://www.3pointproducts.com)**