



3pp[®] ThumSpica[™] Plus



The 3pp[®] ThumSpica[™] Plus is the ideal combination of a soft wrap and a lightweight molded splint that reduces thumb and wrist pain. Used to help treat arthritis, deQuervain's Tenosynovitis, ligament injuries, Gamekeeper's Thumb, sprains and strains.

Features

- The ThumSpica Plus is a soft, wrap-on splint that includes a moldable insert for added control to reduce pain from deQuervain's, arthritis, ligament injuries or after surgery.
- The easy-to-mold insert allows the health care provider to create a custom thumb splint with the comfort of a soft cushioned wrap.
- The molded splint fits into a lined pocket for cushioning and easy application.
- The soft wrap is breathable and hand or machine washable. Latex-free.

Order Information

Measure around the large knuckles of the hand. One size fits most. Product requires fitting by a hand therapist or other medical professional trained in using low temperature thermoplastics.

Product Number	Length
P3019-R (right)	6" - 9"
P3019-L (left)	6" - 9"

Watch our **instructional video on the ThumSpica Plus** on our website: www.3pointproducts.com/product-videos-of-splints/3pp-thumspica-plus-video



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com