



Uno Splint



The Uno splint provides the advantage of a secure wrist immobilizer combined with the convenience of a splint that fits a broad range of hand sizes. Includes stays in the palm and on the top of the hand for full, protective immobilization for tendinitis, wrist sprains, strains and rest for carpal tunnel syndrome, as well as cumulative trauma disorders.

Features

- The innovative two-part design of the Uno fits both the right and left hand and accommodates hand sizes from small to large.
- Once the Uno is adjusted to fit, it is easy to apply and remove. A D-ring strap provides a secure fit around the wrist.
- The Uno is ideal for protecting acutely injured wrists and can be used for children and adolescents.
- Machine washable. Latex-free.

Order Information

One size fits most. 8" in length. Measure circumference at mid-forearm (4" up from the wrist). Fits forearm from 6" to 12". Product Number P3008-U8.

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com