



Wrist Control Splint



Careful attention to detail assures the Wrist Control wrist splint fits properly to relieve the pain and tingling of Carpal Tunnel Syndrome and to rest conditions of the wrist. This splint can also be used for tendinitis, cumulative trauma, wrist sprains and strains.

Features

- Softly cushioned, this splint features a single supporting stay in the palm to protect the wrist and hold it in proper alignment. Use day or night to relieve symptoms of carpal tunnel syndrome.
- The Wrist Control is carefully contoured to avoid pressure on the thumb and allow full finger motion.
- Easy-to-apply D-ring straps provide closure to hold the wrist securely for control and comfort. Three separate straps allow for a custom fit and easy adjustment for changes in wrist circumference throughout the day.
- Machine washable. Latex-free.

Order Information

Measure around the wrist to size.

	Right	Left	Fits
X-Sm	P3004-R1	P3004-R1	4 1/4" to 5 1/4"
Sm	P3004-R2	P3004-R2	5 1/4" to 6 1/4"
Med	P3004-R3	P3004-R3	6 1/4" to 7 1/4"
Lrg	P3004-R4	P3004-R4	7 1/4" to 8 1/4"
X-Lrg	P3004-R5	P3004-R5	8 1/4" to 9 1/4"



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com