

## IMPORTANT SAFEGUARDS & WARNINGS

To ensure that you experience the full range of benefits from the Vanity Planet Buff, we recommend the following:

- Make sure to only use the device on your feet and heels.
- The Smoothing Roller Head is recommended to be used by one person and should not be shared with others.
- Do not use if you have poor circulation or diabetes.
- Before use, consult with your physician if you have any skin conditions or medical issues.
- Do not apply to skin for more than 2-3 seconds at a time.
- Do not use on open wounds, sores, warts, or other problematic areas.
- Overuse may cause skin irritation. Immediately discontinue use if you experience any pain and seek medical attention if skin becomes severely irritated.
- This device should not be used by children or persons with limited physical, sensory, or mental capabilities.
- Store device in a cool, dry place, not in direct sunlight.
- Do not turn on the device if smoothing head is not attached.
- Keep protective cover on the device when not in use.
- Do not disassemble the device.

## ONE-YEAR WARRANTY

This Vanity Planet product is warranted to be free of manufacturer's defects in materials or workmanship for one year from the date of purchase. Damage or wear resulting from an accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair are not covered under this warranty.

Should this product require warranty service (or replacement at our discretion) please contact client service to obtain a Return Merchandise Authorization number (RMA) and return instructions; proof of purchase is required. Products returned without a Vanity Planet generated RMA number will not be accepted and the sender will not receive a refund, replacement, or repaired product.

VANITYPLANET.COM  
INFO@VANITYPLANET.COM  
1.888.870.2132

There are no express warranties except as listed above.  
This warranty gives you specific legal rights and you may have other rights which vary from state to state.

Please do not return this product to the retailer. Contact us directly for warranty service or repair, thank you.

V|P

BUFF

smoothing & refining  
callus remover

#VPBEAUTY | VANITYPLANET.COM


# HEY THERE GOOD LOOKIN'

In this box, you'll find a top notch product waiting patiently for your attention. Let's be honest, you know an awesome buy when you see one, so thanks for choosing our goods to add to your stellar beauty routine.


The way we see it, we're here to help you enhance your already stunning looks. Why? Because we think confidence looks best on everyone.

XOXO  
VANITY PLANET

DON'T FORGET TO SHOW OFF YOUR GOODS BY USING THE  
#MYVP AND #VPBEAUTY HASHTAGS.



- A: PROTECTIVE COVER  
 B: SMOOTHING ROLLER HEAD  
 C: RELEASE BUTTON  
 D: SAFETY SWITCH LOCK  
 E: SLIDE ON/OFF SWITCH  
 F: BATTERY COVER  
 G: BATTERY COVER LOCK  
 H: CLEANING BRUSH



## YOUR BUFF

The Vanity Planet Buff is an innovative pedicure device designed to smooth dry, rough, and callused skin. The Smoothing Roller Head spins 360 degrees at 40 times per second to gently and effectively remove the tough skin on your heels. Faster, safer, and easier to use than traditional scrapers, Buff uses a spinning textured roller rather than metal blades to achieve desired results. This user-friendly device is battery-operated and fits comfortably in your hands, so you can say goodbye to rough patches without the hassle.

For hygienic reasons, the VP Buff is recommended for use by only one person and should not be shared.

## USING YOUR BUFF

Wash and dry your skin before using your Buff. Make sure that the Smoothing Roller Head is placed into the device. Remove the protective cover and prep for use. Turn on the device by simultaneously pressing the safety lock in the center of the On/Off switch and pushing upward. Gently press the Buff against the rough skin on your feet for 2-3 seconds at a time, rotating as needed. Lightly glide the device over the desired areas making sure not to push down with too much force. If pressed too hard, the unit will stop functioning. After use, check to see if the desired degree of softness has been achieved. If not, go over the same area for 2-3 more seconds. Continue to use this process until your skin is silky smooth. Turn off the device when finished. Rinse off the enhanced area or use a damp towel to get rid of excess dead skin. Make sure to clean your Buff after each use. Keep the protective cover on the device when not in use.

NOTE: Only use on your feet and do not use on any broken or bleeding skin.

### CHANGING THE SMOOTHING ROLLER HEAD:

Depending on the frequency of use and size of treated area, the Smoothing Roller Head may need to be replaced every few sessions. Make sure the unit is turned off before replacing the head. Hold the Smoothing Roller Head with one hand and press the release button on the side of the device to eject the head. Place a new Smoothing Roller Head into the device. You will hear a click when it is firmly in place.

## CARING FOR YOUR BUFF

To keep your Vanity Planet Buff in optimum condition, do not use abrasive cleaners or chemicals to clean any part of the device. Clean with water after each use.

### CLEANING YOUR DEVICE:

- Remove the protective cover from the head.
- Clean the device with the included brush; do not switch on when brushing.
- Rinse the unit with water, keeping the head running for 5-10 seconds.
- Turn off the unit and wipe down the body of the device with a clean, dry cloth.
- Replace the protective cover when finished cleaning.

NOTE: The Buff has a water resistant design and can be rinsed under water for cleaning use.

### BATTERY INSTALLATION:

- Remove the batter cover by turning it 90 degrees to the left until it clicks into place, then pull down.
- Insert 2 AA alkaline batteries in the battery compartment.
- Close the battery cover by turning it 90 degrees to the right until the cover locks into place.