



NPWT Dressings:

- **UNI Hydrophobic Dressings**
- **UNI Hydrophilic Dressings**
- **UNI Silver Dressings**

Equinox Medical LLC. specially designed and formulated a hydrophobic environment providing better comfort and less trauma to a patient during dressing changes. Our proprietary dressings accommodate clinician options in treating a wide range of wound types and sizes.

Our foam dressings are made out of high quality reticulated polyurethane/polyethylene foam with an open pore structure (400-600 microns). The use of our hydrophobic, hydrophilic and silver foam dressings were designed with flexible contours for deep or irregularly shaped wounds. UNI foam dressing kits promote cellular activity and formation of granulation tissue in wounds.

The UNI foam dressings help draw wound edges together, remove wound fluids and infectious materials. Our foam technology works actively and can easily adapt to any NPWT pump system.

UNI Black Foam Dressing

- Hydrophobic
- Open-pore, reticulated
- High tensile strength eliminates fraying and reducing contamination and infection

UNI White Foam Dressing

- Hydrophilic
- High density
- Ideal for tunneling and undermining wounds
- Flexible wet or dry, minimizing damage to wound bed

UNI Silver Dressings

- Silverlon antimicrobial contact layer, non-adherent to the wound bed
- Available for 3 day use with NPWT application
- Accelerate wound healing and helps eliminate infection.

***UNI Abdominal now available**



Patient Outcome Driven Technology

UNI Foam Kit Contents

- 1 UNI Foam
- 1 Dome Pad with canister tubing and adapter
- 1 Transparent Drape (1 in Medium Kits, 2 in Large Kits)
- 1 Tubing Clamp



Dome Pad

- Sterile plastic tubing with double pressure lumens that connect the therapy unit to the dressing.
- Double lumen designed to provide consistent NPWT
- Designed to prevent less tube blockage, kinks, pooling and alarm issues
- Dome pad provides a pressure regulating technology at the wound site to provide accurate delivery of prescribed therapy.

Transparent Drape

- Designed with 1, 2 3 step tabs
- Drape can be applied whole or cut to cover the peri-wound area to prevent maceration.

Y-connector

- Y-connectors used for multiple wound applications
- Y-connectors can be adapted to any NPWT pump system

Description	Qty	U/M	Part Number
UNI Hydrophobic Small Kit (10x8x3cm)	20	Case	UN1083S
UNI Hydrophobic Medium Kit (20x13x3 cm)	20	Case	UN20133M
UNI Hydrophobic Large Kit (25x16x3cm)	20	Case	UN25163L
UNI Silver 3 Day Kit Medium (20x13x3 cm)	20	Case	UNSILM
UNI Silver 3 Day Kit Large (25x16x3cm)	20	Case	UNSILL
UNI White Foam Dressing (15x10x1cm)	10	Case	UNIW3001
UNI White Foam Kit w/ Dome and Drape	10	Case	UNIW3000
UNI Adapter	20	Case	UNAD11
Y-Connector	10	Case	YCNN
UNI Dressing Drape Only	10	Case	DRP
UNI Dome Pad/Suction Bell Only	10	Case	USB65

Equinox O2 Medical, LLC.

310 George Washington Hwy, Suite 200, Smithfield, RI 02917

Tel: 800-504-3150 Fax: 800-339-4901 www.EquinoxO2.com

Distributed by:

