

# SILVERLON®

safe, strong, simple™

## Combat Medic Antimicrobial Barrier Wound Dressing

CAT. NO. CM-333-Z

DRESSING SIZE: 3" x 33" (Folded 3" x 4")

PACKAGE DIMENSION: 5 3/4" x 7 3/4" x 1/4"

FDA 510(K) K150256

CAGE CODE: 55U2801

NSN: 6510-01-659-7569

DAPA PRICE: \$56.43



**GENERAL DESCRIPTION:** The Silverlon® Combat Medic Antimicrobial Wound Dressing is a flexible, elastomeric, sterile, silver-plated dressing that delivers antimicrobial silver ions in the dressing when activated by moisture. The silver ions in the dressing kill wound bacteria held in the dressing and provide an antimicrobial barrier for bacterial penetration of the dressing which may help reduce infection.

**MATERIAL COMPOSITION:** The Silverlon® Combat Medic Antimicrobial Wound Dressing is comprised of a single layer of woven 87% nylon and 13% denier elastane that is circumferentially coated with 5.8mg/cm<sup>2</sup> of 99% pure metallic silver and 1% silver oxide.

**INDICATIONS:** Silverlon® Combat Medic Antimicrobial Wound Dressings are indicated for management of a wide variety of wounds including partial and full-thickness wounds, pressure ulcers, venous ulcers, diabetic ulcers, traumatic and surgical wounds, 1st and 2nd degree burns, and donor and graft sites.



Argentum Medical




# Silverlon® Safe. Strong. Simple.

## Simple 2-Dressing Protocol


- Apply Silverlon® pad dressing in operating room
- Apply Silverlon® island dressing at discharge (to be worn for up to 7 days)

### Silverlon® Wound Contact Dressings




PRODUCT #	SIZE	NSA#
WCD-44	4" x 5"	6510-01-522-2168
WCD-1012	10" x 12"	6510-01-509-6925
WCD-466	4" x 66"	6510-01-511-3183

### Silverlon® Burn Contact Dressings



PRODUCT #	SIZE	NSA#
BCD-816	8" x 16"	6510-01-560-5594
BCD-1616	16" x 16"	6510-01-560-5581

### Silverlon® Elastic Burn Wraps



PRODUCT #	SIZE	NSA#
BWD-466	4" x 66"	6510-01-522-2173
BWD-6108	6" x 108"	6510-01-511-3629

### Acute Burn Gloves




PRODUCT #	SIZE	NSA#
ABG-01S	Small	6510-01-560-5965
ABG-01M	Medium	6510-01-522-0305
ABG-01L	Large	6510-01-522-0304

### Silverlon® Digit Sleeve Tapered



PRODUCT #	SIZE	NSA#
DS-112	1.5" to 1" x 12" tapered	6510-01-587-8332

### Combat Medic Wound Care Dressing



PRODUCT #	SIZE	NSA#
CM-333-Z	3" x 33"	6510-01-659-7569

For more details and to read clinical studies on Silverlon Military and EMS Trauma Dressings, visit <http://www.silverlon.com/trauma>

2571 Kaneville Ct.  
Geneva, IL 60134  
888-551-0188  
630-232-2507  
630-232-8005 (FAX)  
[silverlon.com](http://silverlon.com)



**SILVERLON®**  
safe, strong, simple™