

New SCD Pumps!!!

Arjo Flowtron ACS900



The Convenient choice

Flowtron ACS900
Active Compression System

Easy as 1. 2. 3.

APPLY



1. Simply apply
Selected garment(s) to patient

CONNECT



2. Connect garment(s) to pump
Smartsense Auto Garment
Recognition does the rest

START



3. Start Therapy
Simple one button press to start
effective active compression

Cleaning Process:

1. Kept on unit
2. Machines wiped down by EVS
3. Tagged placed
4. Extras stored in Central Sterile Supply

