



Pride in Protection

FINDING THE RIGHT MASK FOR YOU



GHEST **FILTRATION** AND OPTIMAL REATHABILITY. MAKES A LEVEL 3 FEEL AS **COMFORTABLE** AS LEVEL 1.



CELLULOSE **BLEND** OR BI COMPONENT INNER LAYER IS SIMPLY SOFT



DO

NEW MASK AFTER EACH PATIENT





EARLOOPS ARE JUST THE RIGHT LENGTH FOR A PERFECT FIT



MASK WORN COLOUR OUT WITH LOGO ON BOTTOM LEFT FOR PROPER FIT



AS4381:2015 LEVEL





AS4381:2015



AS4381:2015 **REQUIREMENTS**

Bacterial filtration efficiency (BFE)



Sub-micron particulate efficiency at 0.1 µm (PFE) Resistant to penetration by synthetic blood



Differential Pressure (AP)



Flame Spread







@120 mm Hg

< 5.0 mm H₂O/cm²

Class 1

LEVEL

≥ 98% ≥ 98%

@160 mm Hg

< 5.0 mm H₂O/cm²

Class 1



masks are ideal for:



LOW RISK ORAL EXAMS





NO SPRAY ORTHODONTIC TREATMENTS





masks are ideal for:



MEDIUM RISK ORAL EXAMS



MODERATE LEVEL OF SPRAY



WHITENING



HIGH OR LOW-SPEED HANDPIECES



LAB WORK WITH EXPOSURE TO BLOOD DENTAL IMPLANT



masks are ideal for:



HIGH RISK ORAL EXAMS



HIGH LEVEL OF SPRAY



AIR WATER SYRINGE



ORAL SURGERY

GUNZ DENTAL AUSTRALIA

Tel: 02 9935 6677 | www.gunz.com.au | Email: dentalsales@gunz.com.au Unit 4 / 26-34 Dunning Avenue, Rosebery NSW 2018 Australia

