

LabScheduler

Test Request and Scheduling Tool

AND
Discover Precision

Automate test scheduling from initial request to the final results


Time is a precious resource in today's test lab, so the more activities that can be automated, the better. LabScheduler is part of LabCentral, allowing an engineer to submit a test request, which is then managed by lab personnel through fulfillment. Throughout the process, LabScheduler triggers personnel who have an action to perform, and provides feedback to the original test requester.

The requester uses LabCentral to initiate the request, which should include both the parts to be tested and the processes that should be used. These items can either be picked from existing lists, or added at the time the test is being requested.

Once a request is submitted, it enters the scheduling and approval process. The original requester will receive email notifications at each stage in the process so that he can follow the progress through the lab.

An approved test is forwarded to the scheduler to be added to the queue for the appropriate test station. Feedback is available as the test is running via a Test Request display or standard LabCentral dashboards. Once a test is complete, all data and reports are uploaded to the project in LabCentral and accessible via a link in the request or normal LabCentral search methods.

LabScheduler calendar view in LabCentral




Highlights

- Manage, approve, create and view test requests
- Assign tests to workstations remotely via LabCentral
- When combined with standard LabCentral applications, provides a complete test management system for:
 - Users
 - Specimens
 - Test process definitions
 - Test result data
- Keeps test requester in the loop via automated emails throughout the process
- Prevents requests from being queued onto stations without the correct equipment to execute the defined test


LabScheduler

Test Request and Scheduling Tool

Test Request Process Overview



Parts and test processes are selected for each new test request, which is then queued up within LabCentral.



In iTest Battery, tests can be defined and assigned to specific circuits within the test system via the Control Circuits Monitor.

Americas

A&D Technology, Inc.
Ann Arbor, MI USA
Ph: +1 (734)973 1111
www.AandDTech.com

Europe

A&D Europe, GmbH
Griesheim, Germany
PH: +49 (0)6155-605 250
www.AandDEurope.com

A&D Europe - UK Branch
Abingdon, Oxon, UK
PH: +44 (0)1235-550 420
www.AandDEurope.com

Asia

A&D Technology Trading Co.
Shanghai, China
PH: +86 (0)21-3393 2340
www.AandDTech.com

A&D Company, Ltd.
Tokyo, Japan
PH: +81 (0)3-5391 2753
www.AandD.co.jp

Australia

A&D Australasia Pty Ltd.
Kensington, VIC, Australia
PH: +61 (0)3-9372 1522
www.andaustralasia.com.au

