

One-step easy extraction

- Easy to use
home extraction
- Hygienic. Limited
handling of faeces
- No weighing, pipetting or
centrifugation needed
- Stable for 7 days at
ambient temperature
- Tube fits straight into
Dynex DS2 ELISA
automated system



TOTAL SOLUTION for Calprotectin in faeces



CALPRO AS




Arnstein Arnebergsvei 30
N-1366 Lysaker
Norway

Tel: +47 400 04 279
mail@calpro.no
www.calpro.no





CalproSmart Office

- Rapid test for lab use
- Easy access
- Test results in 18min



CE 088



Fully automated Calprotectin ELISA analysis

- Only one dilution needed
- Wide range quantitative ELISA (25 -2500mg/Kg)
- Ideal for analysis of large number of samples (80 samples in 130 min)
- Fully automated system with protocol available for Dynex DS2 ELISA.
- Possibility to store standard curve
- Small and userfriendly kit size
- CalproLab™ ELISA, CE-marked, product No. CALP0170

