


## Certificate of Analysis

### IMMUNOSCAN CCPlus®


REF RA-96 PLUS, RA-96 PLUS RUO

**LOT**            ZS 2177                        2020-10-01

Reagents			LOT		Specification	Result
Ag			ZS 2117	2021-08-09		
CAL	A		ZS 2259	2021-08-26	OD	2.15
CAL	B		ZS 2260	2021-08-26	OD	1.75
CAL	C		ZS 2261	2021-08-26	OD	1.04
CAL	D		ZS 2262	2021-08-26	OD	0.37
CAL	E		ZS 2263	2021-08-26	OD $\geq 0.15$	0.20
CONTROL	REF		ZS 2264	2021-08-26	OD $\geq 0.15$	0.20
CONTROL	+		ZS 2266	2021-08-26	180-340 U/mL	251
					Ratio 4.2-7.6	6.0
CONTROL	-		ZS 2265	2021-08-26	OD	0.02
					Ratio $< 0.95$	0.11
DIL			ZS 2271	2021-08-25	OD	0.05
CONJ			ZS 2267	2021-08-27		
BUF	WASH	20X	ZS 2269	2021-08-25		
SUBS		TMB	ZS 1305	2020-10-01		
H <sub>2</sub> SO <sub>4</sub>		0.5M	ZS 2274	2021-08-25		

It is hereby certified that this LOT has been manufactured, tested and approved according to our quality procedures.

Date of approval: 2019-09-23

  
 \_\_\_\_\_  
 Quality Approved



**Svar Life Science AB**

**Postal address**

P.O. Box 50117  
 SE-202 11 Malmö  
 Sweden

**Visiting address**

Lundavägen 151  
 SE-212 24 Malmö  
 Sweden



[www.svarlifescience.com](http://www.svarlifescience.com)