


New

Xmatrix[®] MINI

eFISHiency Workstation



All-in-One

FISH, *in situ* PCR, ISH

Affordable System

Ten Independent Thermal Cyclers

Apply Micro-reagents for Cost Saving

SMS Prompt

Xmatrix[®] MINI

eFISHiency Workstation




Compact design with
10 slide capacity




Load Slides



Apply hydrophobic barrier
and sealer



Micro-reagent
Dispensing




Place coverslip to make
micro-chamber




View Slides

FISH




- A. FISH EGFR
- B. FISH Break apart ALK
- C. FISH HER2

ISH



- A. HPV Probe
- B. Kappa
- C. Cyclin D1

in situ PCR



- A. HPV 16 *in situ* PCR Ca. Cervix
- B. HPV 16 Negative Control
- C. GAPDH *in situ* PCR Normal Breast

Features & Benefits

- High Performance *in situ* PCR and FISH
- Hybridizer with eXACT™ temperatures
- 10 Independent Thermal Cyclers
- Built-in Touch Screen Display for easy operations
- Compact aesthetics for easy fit in lab
- Facility of on-board wash with effective waste drainage system
- Audio-visual alerts
- GUI with real time display of heating and cooling cycle information
- User-friendly software with ability to add/delete cycles, store protocols for future use and perform up to 45 PCR cycles
- Retrievable run data with date, time and experimental details



Contact Information:

USA: 49026 Milmont Drive, Fremont, CA 94538, USA

India: SEZ Unit, Plot No.4, APIIC SEZ for Aerospace & Precision Engineering
Adibatla, R. R. Dist. - 501510, Telangana State, India

China: Unit 511B, 5F Front Hall, Shanghai Exhibition Center
No.1000 Middle Yan An Road, Shanghai 200041, PRC, China

Italy: Via Savoldini, n.73, 24057, Martinengo, Bergamo, Italy

T: +1-800-421-4149
F: +1-510-824-1490

T: +91-40-27185500
F: +91-40-27185511/21

T: +86-21-2220-6680
F: +86-21-2220-6679

T: +39-036-398-6463



www.biogenex.com

e-mail: customer.service@biogenex.com