eazyplex® SuperBug complete

MDR, XDR, PDR? Carbapenem resistances are increasing!




- Lyophilized tests
- 2 min sample pretreatment
- Portable device
- Real-Time-Detection

Direct screening from swab with eazyplex $^{\circ}$ SuperBug complete A or B - in 30 min! Or confirm your cultures in 20 min!

The test detects the following Carbapenemases as version A:

and now also as version B with OXA-181:

- KPC
- NDM
- OXA-48
- VIM
- OXA-23 group
- OXA-40 group
- OXA-58 group



- KPC
- NDM
- OXA-48
- VIII
- OXA-23 group
- OXA-40 group
- OXA-181

Find out more and contact us under support@amplexbiosystems.de



The **eazyplex**® test system is a freeze dried ready to use amplification reaction and can be stored direct at the workplace. During the amplification the result is shown in Real-Time.

The purchase price of this product includes a non-transferable license under European Patents EP1020534, EP1873260, EP2045337, EP2287338, EP1275715, EP2253717, EP1327679, EP2031058 or their foreign counterparts, owned by Eiken Chemical Co., Ltd.



Amplex BioSystems GmbH

Kerkrader Straße 11 35394 Gießen info@amplexbiosystems.de www.eazyplex.com Bestellservice/Ordering:

Tel.: +49 (0) 80 73/9 16 93-50 Fax: +49 (0) 80 73/9 16 93-55 bestellservice@amplexbiosystems.de Support:

Tel.: +49 (0) 80 73/9 16 93-50 Fax: +49 (0) 80 73/9 16 93-55 info@amplexbiosystems.de

