

Imaging tips

Follow the tips below to take high quality images that are clear, focused, and centered.


OVERVIEW/CLINICAL IMAGES



Take images in a well lit area



Clean your mobile device's lens



Take the photo ~15 cm away from the skin



Tap the mobile screen to focus the image




Take more than one image of your skin condition

DERMOSCOPIC IMAGES



Clean your MoleScope's lens




Ensure the MoleScope fits snugly over your phone's camera. It will either slide down (A) or to the right (B) to align with your camera




Ensure your MoleScope light is ON



Use a fluid media if needed (i.e. a drop of water, oil, gel, alcohol)



MoleScope needs to be in direct contact with skin



Tap the mobile screen to focus the image



Clean MoleScope with an alcohol pad between uses

Note: If you are using another dermatoscope, please follow the instructions provided by the manufacturer.