

# A Simply Smarter Solution for Probe Storage

## Protection for Staff and Patients

New Clean Ultrasound Probe Covers protect intracavity and surface ultrasound probes from the risk of recontamination through handling or the environment after high level disinfection (HLD).

## Custom Designed for Ultrasound Probes

Manufactured in a certified clean room, the covers are custom designed for ultrasound probes. They protect the entire probe including the handle, which should also be disinfected and protected from recontamination as it is routinely handled by staff.


## Correct Storage is a Crucial Reprocessing Step

The Joint Commission expects organizations to use evidence-based national guidelines, such as the ANSI/AAMI ST58:2013 Chemical Sterilization and High-level Disinfection in Health Care Facilities. This standard includes the recommendation that “high level disinfected items should be protected from contamination until the item is delivered to the point of use”.


Clean Ultrasound Probe Covers provide a simple way to deliver best practice patient care while meeting the requirements of reprocessing guidelines and recommendations.



Clean Ultrasound Probe Covers protect the entire probe, including the handle, from recontamination after HLD.



# Point of Care Solution for Fast Workflow



## Reassurance for Your Patients

Patient reassurance is important. When your patient sees you remove a probe from a Clean Ultrasound Probe Cover with a clean HLD label, they can be confident that their exam is being conducted with a probe that has been high level disinfected.

## Clean Ultrasound Probe Covers – The Facts

- Custom designed to fit most intracavity probes and many surface probes
- Protect the entire probe including the handle up to the strain relief
- Enclosed wire wrap fastens cover around the probe
- Internal cover size 4.3" W x 13.6" L
- Manufactured in a certified clean room (ISO 14644-1, Class 8)
- Latex free
- Single use
- Available in boxes of 100

### CDC and FDA Health Alert<sup>1</sup>

The Centers for Disease Control and Prevention (CDC) and the FDA are now alerting healthcare providers and facilities about the need to properly maintain, clean, and disinfect or sterilize reusable medical devices.

All healthcare facilities with reusable medical devices are being urged to review their practices to ensure they comply with device manufacturers' instructions and current standards and guidelines.

<sup>1</sup> CDC Health Alert Network September 11, 2015. Available from: <http://emergency.cdc.gov/han/han00382.asp>

Call **1-844-TROPHON**

and quote product code N00102 to order



#### Nanosonics Limited (Manufacturer)

14 Mars Road  
Lane Cove NSW 2066  
Australia  
T +61 2 8063 1600  
E [info@nanosonics.com.au](mailto:info@nanosonics.com.au)  
[www.nanosonics.com.au](http://www.nanosonics.com.au)

#### Nanosonics, Inc. (Distributor)

11793 Technology Lane  
Fishers, IN 46038  
USA  
T 1-844-TROPHON 1-844-876-7466  
E [info@trophon.com](mailto:info@trophon.com)  
[www.trophon.com](http://www.trophon.com)

#### Nanosonics Europe Limited (Distributor)

Unit 7, Linfit Court, Colliers Way  
Clayton West, Huddersfield HD8 9WL  
United Kingdom  
T +44 1484 860581  
E [ukinfo@nanosonics.eu](mailto:ukinfo@nanosonics.eu)  
[www.nanosonics.eu](http://www.nanosonics.eu)

#### Nanosonics Europe GmbH (EU Representative)

Poppenbueteler Bogen 66  
22399 Hamburg  
Germany  
T +49 40 46856885  
E [info@nanosonics.eu](mailto:info@nanosonics.eu)  
[www.nanosonics.eu](http://www.nanosonics.eu)