

trophon®

Companion Cleaning Wipes

Cleaner and disinfectant

Multi-purpose cleaning and low level disinfection (LLD) wipes, uniquely designed for ultrasound transducers and equipment at point of care.

Applications:




Monitors, touchscreens and keyboards



Non-critical ultrasound probes



Probe holders



Plugs



Cables




trophon multi-purpose cleaning and disinfection wipes replace detergent products. The active ingredients effectively clean and deliver low level disinfection efficacy in accordance with international standards.

nanosonics****
Infection Prevention. For Life.



Ready to use – no need to add chemicals or mess around with pump canisters



Fast and efficient disinfection of ultrasound equipment including the console, monitors, keyboards, cables and probe holders

Effective against viruses, bacteria and fungi including*:

- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Staphylococcus aureus*
- *Acinetobacter baumannii*
- *Enterococcus faecalis* – vancomycin resistant (VRE)
- *Klebsiella pneumoniae* (NDM-1)
- *Staphylococcus aureus* – methicillin-resistant (MRSA)
- *Candida albicans*
- Hepatitis B virus
- Hepatitis C virus
- Herpes simplex virus type 1
- Herpes simplex virus type 2
- HIV-1
- Influenza A virus H1N1
- Norwalk virus (norovirus)
- Rotavirus
- SARS associated coronavirus (SARS)
- Vaccinia virus

*This is just a subset of organisms that are inactivated by trophon Companion Wipes. Please refer to product information.



trophon Companion Cleaning Wipes are also the perfect partner for one step cleaning to prepare ultrasound probes for high level disinfection (HLD).



trophon Companion Drying Wipes are highly absorbent to complement the Cleaning Wipes and are used to dry probes and surfaces before and after LLD and HLD, leaving surfaces completely dry.

A partner you can trust

Nanosonics is a global innovator in infection prevention. trophon2, the latest innovation in ultrasound probe HLD, features an enhanced design; simple and fast workflows; and all new AcuTrace™ for digital record keeping and seamless integration with your hospital IT system.

Description	Code	Qty
trophon Companion Cleaning Wipes	N00112-US	6 x 100 wipes
trophon Companion Drying Wipes	N00113	6 x 100 wipes
trophon Companion Wipes Wall Bracket	N00114	1

To learn more visit www.nanosonics.us or call 1-844-TROPHON (1-844-876-7466)

Nanosonics Limited (Manufacturer), 14 Mars Road, Lane Cove NSW 2066, Australia MM0811-US-BR V05