

# Capp Harmony

## SEROLOGICAL PIPETTES

Capp Harmony pipettes are made from transparent medical-grade polystyrene sterilized by gamma irradiation. Each lot is tested for endotoxins (C.5 Eu/ML) (ISO TO993-11) and Cytotoxicity (ISO 10993-5) and are guaranteed non-haemolytic (ISO 10993-4). Ideal for use within tissue culture applications.

### Features:

- Designed for accurate and rapid fluid dispensing, these polystyrene pipettes are sterile, pyrogen-free and cotton-plugged for a variety of liquid handling applications
- The sharp, black, highly legible permanent graduations prevent ambiguity and guarantee precise reading and accurate dispensing to +1-2%.
- Pipettes are supplied individually wrapped in fibre-free colour-coded paper and plastic combination that is easy to open
- Each pipette comes with an effective filter barrier to prevent overflow of liquid into pipette
- Reverse graduations show the volume removed or remaining on all pipettes larger than 1mL.
- Generous negative graduations provide additional pipetting volume.
- The uniform mouth piece fits all major pipettors and is colour-coded for ease of identification

## Ordering Information

Cat. No.	Description	Volume	Package	Colour
SP-1-C	Serological pipette	1mL	500 pcs.	Yellow
SP-2-C	Serological pipette	2mL	500 pcs.	Green
SP-5-C	Serological pipette	5mL	200 pcs.	Blue
SP-10-C	Serological pipette	10mL	200 pcs.	Orange
SP-25-C	Serological pipette	25mL	150 pcs.	Red
SP-50-C	Serological pipette	50mL	100 pcs	Purple
SP-100-C	Serological pipette	100mL	50 pcs.	Pink

