


CappWash ELISA Plate Washer Features



Easy to Operate

Hold the manifold towards the plate and press the button to dispense the washing solution. Release the button to aspirate the liquid, which will go directly to the waste collection bottle.

Simple & Robust Design

The 2 in 1 tube design prevents the channels from clogging. Hand held work allows the user to perform a more qualified washing of the plate, compared to automated washers.

Manual ELISA Plate Washer

Ideal tool for low scale and reliable ELISA plate washing. Simple to use and economical, users can choose to simply order the manifold or the complete set, which includes a vacuum pump, bottles and tubing.