



Rotor A-2-MTP

- Max. speed: 2,204 × g (4,680 rpm)
- 2-place swing-bucket rotor for MTP, PCR and Deepwell Plates
- Catalog No. 022654403



Rotor F-45-18-17-Cryo

- Max. speed: 8,324 × g (8,900 rpm)
- 18-place rotor for Cryovials®
- Catalog No. 022654161



Rotor FA-45-24-11-HS

- Max. speed: 30,130 × g (17,500 rpm)
- Max. capacity: 24 × 1.5/2.0 mL tubes
- Aerosol-tight
- PTFE coated for increased chemical resistance
- Catalog No. 022654080



Rotor F-45-30-11

- Max. speed: 20,817 × g (14,000 rpm)
- Max. capacity: 30 × 1.5/2.0 mL tubes
- PTFE coated for increased chemical resistance
- Catalog No. 022654004



Rotor F-45-64-5-PCR

- Max. speed: 13,543 × g (11,800 rpm)
- PCR strip rotor for up to 64 × 0.2 mL PCR tubes
- Catalog No. 022654209



Rotor F-35-6-30

- Max. speed: 7,745 × g (7,830 rpm)
- Rotor for 6 × 15/50 mL conical, 12 × Vacutainer® tubes
- Catalog No. 022654306



Rotor FA-45-16-17

- Max. speed: 21,191 × g (14,200 rpm)
- Max. capacity: 16 × 5.0 mL Tubes
- Aerosol-tight
- Catalog No. 5427750002



Rotor FA-45-24-11-Kit

- Max. speed: 19,090 × g (13,200 rpm)
- Max. capacity: 24 x spin columns
- Aerosol-tight
- Catalog No. 022654128



Rotor FA-45-30-11

- Max. speed: 20,817 × g (14,000 rpm)
- Max. capacity: 30 × 1.5/2.0 mL tubes
- Aerosol-tight
- PTFE coated for increased chemical resistance
- Catalog No. 022654047



Rotor FA-45-48-11

- Max. speed: 18,210 × g (12,700 rpm)
- Max. capacity: 48 × 1.5/2.0 mL tubes
- Aerosol-tight
- Catalog No. 5427754008



Rotor F-45-48-11

- Max. speed: 18,210 × g (12,700 rpm)
- Max. capacity: 48 × 1.5/2.0 mL tubes
- Catalog No. 5427755004



Rotor S-24-11-AT

- Max. speed: 16,049 × g (12,700 rpm)
- Swing-bucket rotor for 24 × 1.5/2.0 mL tubes
- Aerosol-tight
- Catalog No. 5427757007