



COBRA[®]

Pipette



VWR  | vwr.com
We Enable Science

Biotix[®]
Fluid Innovation

DESIGNED FOR COMFORT AND ACCURACY.

The Cobra Manual Pipettor System combines air displacement technology with the latest ergonomic research for the most comfortable pipetting imaginable. The low pipetting forces required reduce the risk of repetitive stress disorder, and the built-in control settings ensure accuracy in ranges from 0.2 μL to 1000 μL .

ERGONOMICS

- 1 Easily used by left-hand or right-hand users thanks to the adjustable tip-ejector button
- 2 Extra-long finger hook rests securely into hands of all sizes for a firm grip

CONSTRUCTION

Light-weight design falls within the ideal weight range

- 3 Stainless-steel piston provides precision and durability
- 4 Fully autoclavable at 121° C, 1 atm for 20 minutes

Full 12-month warranty

PERFECT FIT

- 5 All Biotix™ Manual Pipette Tips feature FlexFit™ for the best fit with the least amount of force.

FlexFit technology provides flexibility on the proximal end of tips, reducing the necessary insertion and ejection forces between pipettor and pipette tips.



Cobra Variable Volume Pipettor Line

VWR CAT. NO.	PRODUCT TITLE	VOLUME RANGE (µL)	SYSTEMATIC ERROR/ACCURACY (%)	RANDOM ERROR/PRECISSION (%)
89525-042	Cobra® Single Channel Pipettor, Biotix	0.2 to 2 µl	±13.2% - ±1.7%	≤6.6% - 0.8%
89525-044	Cobra® Single Channel Pipettor, Biotix	1 to 10 µl	±2.8% - ±1.1%	≤1.3%- 0.4%
89525-046	Cobra® Single Channel Pipettor, Biotix	2 to 20 µl	±5.5%-±1.0%	≤1.7%- 0.3%
89525-048	Cobra® Single Channel Pipettor, Biotix	10 to 100 µl	±3.9%-±0.8	≤1.1%- 0.2%
89525-050	Cobra® Single Channel Pipettor, Biotix	20 to 200 µl	±2.8%-±0.8%	≤1.1%- 0.2%
89525-052	Cobra® Single Channel Pipettor, Biotix	100 to 1000 µl	±3.3%-±0.8%	≤0.7%-0.2%
89525-054	Cobra® 8 Channel Pipettor, Biotix	0.5 to 10 µl	±17.6%-±2.2%	≤8.8%-1.1%
89525-056	Cobra® 8 Channel Pipettor, Biotix	2 to 20 µl	±5.5%-±2.0%	≤4.4%- 0.8%
89525-058	Cobra® 8 Channel Pipettor, Biotix	20 to 200 µl	±2.8%- ±1.1%	≤1.4%- 0.3%
89525-060	Cobra® 8 Channel Pipettor, Biotix	20 to 300 µl	±3.7%- ±1.1%	≤1.3%- 0.4%
89525-062	Cobra® 12 Channel Pipettor, Biotix	0.5 to 10 µl	±17.6%-±2.2%	≤8.8%-1.1%
89525-064	Cobra® 12 Channel Pipettor, Biotix	2 to 20 µl	±5.5%-±2.0%	≤4.4%-0.8%
89525-066	Cobra® 12 Channel Pipettor, Biotix	20 to 200 ul	±2.8%- ±1.1%	≤1.4%-0.3%
89525-068	Cobra® 12 Channel Pipettor, Biotix	20 to 300 ul	±3.7%- ±1.1%	≤1.3%- 0.4%
89525-070	Carousel Stand for 7 pipettes			



©2015, Biotix, Inc. All rights reserved

Biotix, Fluid Innovation, Delta Filter, X-Resin, Chem-Resin, Blade, FlexFit, Statix, Leave Nothing Behind, BioReady, Tip Eject, and XTreme Quality are registered trademarks or trademarks of Biotix, Inc. All other brands and names contained herein are the property of their respective owners.



1.800.932.5000 | vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2015 VWR International, LLC. All rights reserved.