

OXFORD
LAB PRODUCTS

BenchMate® Centrifuge Tubes




Available at

Pipette.com

info@pipette.com • 800-243-3232

Oxford Lab Products
BenchMate® Centrifuge Tubes
Product Overview



Free of Contaminants

Produced in a fully automatic, human touch-free clean room facility ensuring the tubes are free from DNase, RNase, human DNA, endotoxins, pyrogens, PCR inhibitors, lubricants, and heavy metals

Oxford Lab Products

BenchMate® Centrifuge Tubes

Features and Benefits



50 mL Centrifuge Tubes BenchMate® Centrifuge Tubes

Maximum Speed: 17,000 x g
Temperature Range - 80 °C to +121°C

Catalog No.	Description	Quantity
OCT-50B	BenchMate® 50 mL Sterile Centrifuge Tubes with Attached Flat Cap, Bulk Pack	25 Tubes/Resealable Pack, 20 Packs per Case
OCT-50R	BenchMate® 50 mL Sterile Centrifuge Tubes with Attached Flat Cap, Racked	25 Tubes/Rack, 20 Packs per Case



15 mL Centrifuge Tubes BenchMate® Centrifuge Tubes

Maximum Speed: 17,000 x g
Temperature Range - 80 °C to +121°C

Catalog No.	Description	Quantity
OCT-15B	BenchMate® 15 mL Sterile Centrifuge Tubes with Attached Flat Cap, Bulk Pack	25 Tubes/Resealable Pack, 20 Packs per Case
OCT-15R	BenchMate® 15 mL Sterile Centrifuge Tubes with Attached Flat Cap, Racked	25 Tubes/Rack, 20 Packs per Case

Oxford Lab Products
BenchMate®
Additional Products

Original BenchMate®
Pipette



BenchMate® L Single
Channel Pipette



BenchMate® E
Single and Multi Channel Electronic Pipettes



Digital Vortex Mixer
(VM-D)



Microcentrifuge
(C8)



Digital Plate Shaker
(S4P-D)



High Speed Microcentrifuge
(C12V)