


# XA 4Y Analytical Balances


Analytical weighing of the highest quality and accuracy using professional balances



XA 4Y d = 0,01 mg



XA 4Y, d = 0,1 mg



Automatic control of the level



Intuitive operation and touch screen

## Functions

- |                |                       |                             |                                |                   |
|----------------|-----------------------|-----------------------------|--------------------------------|-------------------|
| Parts counting | Percent weighing      | Pipettes calibration        | Cooperation with titrators     | Moveable range    |
| Dosing         | Statistics            | Statistical quality control | Ambient conditions measurement | Replaceable unit  |
| Checkweighing  | Animal weighing       | Autotest                    | GLP procedures                 | Multilingual menu |
| Formulations   | Differential weighing | Density determination       | Proximity sensors              |                   |

## Features

### Reliable Results and Excellent Weighing Parameters

The best weighing parameters enable applying XA 4Y analytical balances in majority of demanding laboratory processes.

### Innovative 2-Point Adjustment

The unique 2-point adjustment system guarantees the best accuracy of the weighing results and minimizes linearity errors thus ensuring reliable measurements within the whole measuring range.

### Significantly Fast Measurements

Powerful processor offers new possibilities of operation assuring short indication stabilization time at a respective repeatability.

### Second to None Repeatability and Compliance with USP

XA 4Y analytical balances feature the highest measurements accuracy, excellent repeatability and are compliant with USP requirements (Chapter 41 and 1251).

### Intuitive Operation and Touch Screen

5.7" colour touch screen enables intuitive operation and easy access to numerous applications and functions of the weighing instrument.

### Automatic Control of the Level

Levelling system facilitates adjustment of device level, it also uninterruptedly controls the level state, and informs about potential level deviations.

### Spacious Weighing Chamber

Large weighing chamber enables convenient operation using laboratory vessels of different dimensions.

### Numerous Options of Data Management

Extensive storage capacity allows record of all measurement data in a form of complex reports and statistical graphs.

# Technical Specifications

	XA 52.4Y	XA 82/220.4Y	XA 110.4Y	XA 120/250.4Y
<b>Maximum capacity [Max]</b>	52 g	82/220 g	110 g	120 g /250 g
<b>Minimum load</b>	1 mg	1 mg	1 mg	1 mg
<b>Readability [d]</b>	0.01 mg	0.01 mg / 0.1 mg	0.01 mg	0.01g /0.1 mg
<b>Verification scale interval [e]</b>	1 mg	1 mg	1 mg	1 mg
<b>Tare range</b>	-52 g	-220 g	-110 g	-250 g
<b>Repeatability(5% Max)*</b>	0.006 mg (Rt ≤ 3g)	0.007 mg (Rt ≤ 10 g)	0.007 mg (Rt ≤ 5 g)	0.007 mg (Rt ≤ 10 g)
<b>Repeatability(Max)*</b>	0.01 mg	0.06 mg	0.02 mg	0.06 mg
<b>Linearity</b>	±0.03 mg	±0.06 / ±0.2 mg	±0.06 mg	±0.06 / ±0.2 mg
<b>Eccentric load deviation</b>	0.03 mg	0.2 mg	0.06 mg	0.2 mg
<b>Sensitivity temperature drift**</b>	$1 \times 10^{-6} / ^\circ\text{C} \times \text{Rt}$	$1 \times 10^{-6} / ^\circ\text{C} \times \text{Rt}$	$1 \times 10^{-6} / ^\circ\text{C} \times \text{Rt}$	$1 \times 10^{-6} / ^\circ\text{C} \times \text{Rt}$
<b>Sensitivity time drift</b>	$1 \times 10^{-6} / \text{Year} \times \text{Rt}$	$1 \times 10^{-6} / \text{Year} \times \text{Rt}$	$1 \times 10^{-6} / \text{Year} \times \text{Rt}$	$1 \times 10^{-6} / \text{Year} \times \text{Rt}$
<b>Minimum weight (U=1%, k=2)</b>	1.2 mg	1.4 mg	1.4 mg	1.4 mg
<b>Minimum weight (USP)</b>	12 mg	14 mg	14 mg	14 mg
<b>Stabilization time</b>	~ 4 s	~ 4 s	~ 4 s	~ 4 s
<b>Adjustment</b>	internal	internal	internal	internal
<b>Moveable range</b>	—	Yes	—	Yes
<b>Verification</b>	Yes	Yes	Yes	Yes
<b>OIML Class</b>	I	I	I	I
<b>Indicator fastening</b>	35 cm cable, wireless connection (option)***	35 cm cable, wireless connection (option)***	35 cm cable, wireless connection (option)***	35 cm cable, wireless connection (option)***
<b>Display</b>	5.7" colour, resistive touch screen	5.7" colour, resistive touch screen	5.7" colour, resistive touch screen	5.7" colour, resistive touch screen
<b>Keypad</b>	8 keys	8 keys	8 keys	8 keys
<b>Protection class</b>	IP 43	IP 43	IP 43	IP 43
<b>Databases</b>	19	19	19	19
<b>Touch-free operation</b>	2 programmable proximity sensors	2 programmable proximity sensors	2 programmable proximity sensors	2 programmable proximity sensors
<b>USB-A</b>	2	2	2	2
<b>Ethernet</b>	10 / 100 Mbit	10 / 100 Mbit	10 / 100 Mbit	10 / 100 Mbit
<b>RS 232</b>	2	2	2	2
<b>Wireless connection</b>	802.11 b/g/n	802.11 b/g/n	802.11 b/g/n	802.11 b/g/n
<b>IN/OUT</b>	4 × IN, 4 × OUT	4 × IN, 4 × OUT	4 × IN, 4 × OUT	4 × IN, 4 × OUT
<b>Power supply</b>	13.5 ÷ 16 V DC	13.5 ÷ 16 V DC	13.5 ÷ 16 V DC	13.5 ÷ 16 V DC
<b>Power consumption</b>	10 W	10 W	10 W	10 W
<b>Operating temperature</b>	+10 ÷ +40 °C	+10 ÷ +40 °C	+10 ÷ +40 °C	+10 ÷ +40 °C
<b>Atmospheric humidity****</b>	40 ÷ 80%	40 ÷ 80%	40 ÷ 80%	40 ÷ 80%
<b>Transport and storage temperature</b>	-20 ÷ +50 °C	-20 ÷ +50 °C	-20 ÷ +50 °C	-20 ÷ +50 °C
<b>Weighing pan dimensions</b>	ø 90 mm open-work ø 85 mm standard (option)*****	ø 90 mm open-work ø 85 mm standard (option)*****	ø 90 mm open-work ø 85 mm standard (option)*****	ø 90 mm open-work ø 85 mm standard (option)*****
<b>Weighing chamber dimensions</b>	168 × 160 × 228 mm	168 × 160 × 228 mm	168 × 160 × 228 mm	168 × 160 × 228 mm
<b>Weighing device dimensions</b>	542 × 206 × 303 mm	542 × 206 × 303 mm	542 × 206 × 303 mm	542 × 206 × 303 mm
<b>Net weight</b>	9.8 kg	9.8 kg	9.8 kg	9.8 kg
<b>Gross weight</b>	14.3 kg	14.3 kg	14.3 kg	14.3 kg
<b>Packaging dimensions</b>	720 × 385 × 485 mm	720 × 385 × 485 mm	720 × 385 × 485 mm	720 × 385 × 485 mm

Rt net weight  
 \* repeatability is expressed as a standard deviation from 10 weighing cycles  
 \*\* parameter determined in the following temperature range: +15 ÷ +35 °C  
 \*\*\* optional solution on purchase order  
 \*\*\*\* non-condensing conditions  
 \*\*\*\*\* ø 85 mm standard weighing pan on purchase order

Values of parameters provided in Technical Specifications table have been determined under stable laboratory conditions. Due to ambient conditions impact or/and balance setup, the above parameters may vary for environments other than laboratory.

	XA 210.4Y	XA 220.4Y	XA 310.4Y
Maximum capacity [Max]	210 g	220 g	310 g
Minimum load	1 mg	10 mg	10 mg
Readability [d]	0.01 mg	0.1 mg	0.1 mg
Verification scale interval [e]	1 mg	1 mg	1 mg
Tare range	-210 g	-220 g	-310 g
Repeatability(5% Max)*	0.007 mg (Rt ≤ 10g)	0.06 mg (Rt ≤ 10g)	0.07 mg (Rt ≤ 15g)
Repeatability(Max)*	0.03 mg	0.1 mg	0.15 mg
Linearity	±0.1 mg	±0.2 mg	±0.3 mg
Eccentric load deviation	0.1 mg	0.2 mg	0.3 mg
Sensitivity temperature drift**	$1 \times 10^{-6} / ^\circ\text{C} \times \text{Rt}$	$1 \times 10^{-6} / ^\circ\text{C} \times \text{Rt}$	$1 \times 10^{-6} / ^\circ\text{C} \times \text{Rt}$
Sensitivity time drift	$1 \times 10^{-6} / \text{Year} \times \text{Rt}$	$1 \times 10^{-6} / \text{Year} \times \text{Rt}$	$1 \times 10^{-6} / \text{Year} \times \text{Rt}$
Minimum weight (U=1%, k=2)	1,4 mg	12 mg	14 mg
Minimum weight (USP)	14 mg	120 mg	140 mg
Stabilization time	~ 4 s	~ 2.5 s	~ 2.5 s
Adjustment	internal	internal	internal
Moveable range	—	—	—
Verification	Yes	Yes	Yes
OIML Class	I	I	I
Indicator fastening	35 cm cable. wireless connection (option)***	35 cm cable. wireless connection (option)***	35 cm cable. wireless connection (option)***
Display	5.7" colour, resistive touch screen	5.7" colour, resistive touch screen	5.7" colour, resistive touch screen
Keypad	8 keys	8 keys	8 keys
Protection class	IP 43	IP 43	IP 43
Databases	19	19	19
Touch-free operation	2 programmable proximity sensors	2 programmable proximity sensors	2 programmable proximity sensors
USB-A	2	2	2
Ethernet	10 / 100 Mbit	10 / 100 Mbit	10 / 100 Mbit
RS 232	2	2	2
Wireless connection	802.11 b/g/n	802.11 b/g/n	802.11 b/g/n
IN/OUT	4 × IN. 4 × OUT	4 × IN. 4 × OUT	4 × IN. 4 × OUT
Power supply	13.5 ÷ 16 V DC	13.5 ÷ 16 V DC	13.5 ÷ 16 V DC
Power consumption	10 W	10 W	10 W
Operating temperature	+10 ÷ +40 °C	+10 ÷ +40 °C	+10 ÷ +40 °C
Atmospheric humidity****	40 ÷ 80%	40 ÷ 80%	40 ÷ 80%
Transport and storage temperature	-20 ÷ +50 °C	-20 ÷ +50 °C	-20 ÷ +50 °C
Weighing pan dimensions	ø 90 mm open-work ø 85 mm standard (option)*****	ø 100 mm standard	ø 100 mm standard
Weighing chamber dimensions	168 × 160 × 228 mm	168 × 160 × 228 mm	168 × 160 × 228 mm
Weighing device dimensions	542 × 206 × 303 mm	542 × 206 × 303 mm	542 × 206 × 303 mm
Net weight	9.8 kg	9.8 kg	9.8 kg
Gross weight	14.3 kg	14.3 kg	14.3 kg
Packaging dimensions	720 × 385 × 485 mm	720 × 385 × 485 mm	720 × 385 × 485 mm

Rt net weight

\* repeatability is expressed as a standard deviation from 10 weighing cycles

\*\* parameter determined in the following temperature range: +15 ÷ +35 °C


\*\*\* optional solution on purchase order

\*\*\*\* non-condensing conditions




\*\*\*\*\* ø 85 mm standard weighing pan on purchase order

Values of parameters provided in Technical Specifications table have been determined under stable laboratory conditions. Due to ambient conditions impact or/and balance setup, the above parameters may vary for environments other than laboratory.



## Dimensions



XA 4Y, d = 0.01 mg



XA 4Y, d = 0.1 mg



## Accessories

### Weighing Tables

- granite antivibration table
- antivibration tables for laboratory balances
- professional weighing table
- workstation for pipettes calibration

### Professional Weighing

- Adapter for calibration of XA17 i XA100 series pipettes
- KIT 85 density determination kit

### Ambient Conditions

- DJ-04 anti-static ioniser
- THB-Y ambient conditions module

### Peripheral Devices

- Epson dot matrix printer
- barcode scanners
- WD-5/3Y LCD display

### Cables, Converters

- P0108: RS 232 cable (balance-computer)
- P0167: RS 232 cable (balance-computer)
- P0151: RS 232 cable (balance - Epson printer)

### Electrical Accessories

- ZR-02 power supply with battery

## Dedicated Software

### Label Editor R02

- designing label templates
- sending graphics and fonts to label printers
- printing label templates using connected printers

### R-LAB

- collecting measurements
- carrying out statistical analysis of measurements
- customized graphs and reports

### RADWAG Remote Desktop

- remote operation via computer, mobile phone or tablet
- sending text messages
- version for Windows 10 and Android systems

### Pipettes

- determining measurement errors of pipettes volume
- accordance with ISO 8655
- calibration of single-channel and multi-channel pipettes
- calibration of fixed-volume and variable-volume pipettes