



## Buy 2 cases of CryoELITE® Cryogenic Vials Receive a FREE case of Keep-IT® Freezer Boxes

### CryoELITE® Cryogenic Vials

- Lot certified RNase / DNase and endotoxin free and non-pyrogenic providing assurance for sample integrity
- Unrivaled cap seal exceeds DOT and IATA regulations ensuring ultimate protection of samples during transportation and demanding freeze-thaw handling
- Made from low binding, cryogenic grade virgin polypropylene
- Screw cap can be easily removed with one hand

#### Freestanding


- Loctagon™ Vial Skirt provides stability in freestanding position and locks into CryoELITE® Benchmate Rack in order to provide easy open and close with one hand
- Colored caps allow for color coding projects along with WHEATON colored freezer and storage boxes
- Optional 2D Data Matrix Bar Code Insert provides unique identifier for traceability
- Bottom format allows unrestricted view of 2D bar code for convenient automated scanning

Cat No.	Size (mL)	Color	Writing Patch	Sterile	Dimensions Dia. x H (mm)	Qty / Case
CryoELITE® Vials, Freestanding, Internal Thread						
W985922	2	Natural	Yes	Yes	12 x 50	500
W985916	2	White	Yes	Yes	12 x 50	500
W985917	2	Red	Yes	Yes	12 x 50	500
W985918	2	Pink	Yes	Yes	12 x 50	500
W985919	2	Yellow	Yes	Yes	12 x 50	500
W985920	2	Green	Yes	Yes	12 x 50	500
W985921	2	Blue	Yes	Yes	12 x 50	500
CryoELITE® Vials, Freestanding, External Thread						
W985852	2	Natural	Yes	Yes	12 x 37	500
W985863	2	White	Yes	Yes	12 x 49	500
W985864	2	Red	Yes	Yes	12 x 49	500
W985865	2	Pink	Yes	Yes	12 x 49	500
W985866	2	Yellow	Yes	Yes	12 x 49	500
W985867	2	Green	Yes	Yes	12 x 49	500
W985868	2	Blue	Yes	Yes	12 x 49	500

### Keep-IT® Freezer boxes

- KeepIT® Freezer Boxes provide an ideal method for batching and storing samples
- Six different colors match the colors of CryoELITE® Cryogenic Vials, creating alternatives for batching and identifying groups of samples
- KeepIT®-25 accommodates up to 25 internal or external threaded cryogenic Vial Sizes: 1.2 - 2mL
- KeepIT®-81 accommodates 81 external threaded cryogenic vials (sizes: 1.2 - 2mL)
- Openings in bottom facilitate scanning CryoELITE® 2D Data Matrix Bar Code Inserts
- Made from Eastman Tritan™ BPA free, shatter resistant resin
- Standard footprint compatible with liquid nitrogen storage shelves and freezer drawers

Cat. No.	Color	Dimensions (L x W x H) (mm)	Qty / Case
KeepIT®-25 For External Thread Vials			
W651702-W	White	75 x 75 x 52	10
W651702-R	Red	75 x 75 x 52	10
W651702-P	Pink	75 x 75 x 52	10
W651702-Y	Yellow	75 x 75 x 52	10
W651702-G	Green	75 x 75 x 52	10
W651702-B	Blue	75 x 75 x 52	10
Low Profile KeepIT®-81 For External Thread Vials			
W651703-W	White	130 x 130 x 52	10
W651703-R	Red	130 x 130 x 52	10
W651703-P	Pink	130 x 130 x 52	10
W651703-Y	Yellow	130 x 130 x 52	10
W651703-G	Green	130 x 130 x 52	10
W651703-B	Blue	130 x 130 x 52	10



Available at

**Pipette**.com

info@pipette.com • 800-243-3232

To claim your redemption (redeem by 7/31/19)

1. Obtain your Proof of Purchase (dated between 1/1/19 - 6/30/19). Packing slip, purchase order and invoice are all considered valid for proof of purchase.
2. Visit [www.dwkdeals.com](http://www.dwkdeals.com)
3. Attach your Proof of Purchase and fill-out the online form.
4. Enter the part numbers for freeware. Free items must be of equal or lesser value. Visit [www.DWK.com](http://www.DWK.com) for more details. Not for resale. Limited to two redemptions per promotion, per customer and/or per laboratory.

Offer valid for US and Canadian customers only. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. Offer valid while supplies last. Not valid on prior purchases. DWK distributors and employees are prohibited from participating. Void in government accounts or where prohibited by law or company policy. Prices subject to change without notice. Trademarks contained herein are the property of their respective owners.



DWK Life Sciences  
1501 North 10th Street  
Millville, NJ 08332  
United States

US & Canada: 800.225.1437  
Int'l: 856.825.1100  
Fax: 856.825.1368  
[www.dwk.com](http://www.dwk.com)