



**Buy 2 KIMBLE® Glassware Starter Packs
 Receive a FREE GL 45 Media Bottle Starter Pack**

Low Form Beaker Starter Packs

- Double capacity scale
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

Cat. Number	Capacity (mL)	Type	Case Qty
14080-01	1 each of 50, 100, 250, 600, 1000	Standard Wall	1 Set of 5
14085-01	1 each of 250, 400, 600, 1000	Heavy Duty	1 Set of 4

Erlenmeyer Starter Pack

- The pack consists of 5 flasks in sizes of 50, 125, 250, 500 and 1000mL
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

Cat. Number	Capacity (mL)	Rubber Stopper Size	Case Qty
26520-1	50, 125, 250, 500, 1000	1, 5, 6, 7, 9	1

Class B Cylinders Starter Pack

- "TD" indicates the cylinder is calibrated to deliver
- Class B
- Single metric scale, with bumper
- The pack consists of 5 cylinders in sizes of 10, 25, 50, 100 and 250mL
- Designed from ASTM Specification E1272, Style I, Class B requirements
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements


Cat. Number	Capacity (mL)	Case Qty
20024-01	10, 25, 50, 100, 250	1

GL 45 Media Bottle Starter Pack

- Enhanced graduations and marking spot made with chemically resistant white enamel paint
- 30 mm ID opening
- Autoclavable
- Supplied with blue polypropylene linerless GL 45 screw thread caps
- Replaceable clear (natural) drip-free polypropylene pour ring is included on each bottle
- The pack consists of two 100 mL, three 250 mL, three 500 mL, and two 1000 mL bottles
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

Cat. Number	Capacity (mL)	Case Qty
14395-01	100 (3), 250 (3), 500 (3), 1000 (2)	1





Available at

Pipette.com

info@pipette.com • 800-243-3232

To claim your redemption (redeem by 7/31/19)

1. Obtain your Proof of Purchase (dated between 1/1/19 - 6/30/19). Packing slip, purchase order and invoice are all considered valid for proof of purchase.
2. Visit www.dwkdeals.com
3. Attach your Proof of Purchase and fill-out the online form.
4. Enter the part numbers for freeware. Free items must be of equal or lesser value. Visit www.DWK.com for more details. Not for resale. Limited to two redemptions per promotion, per customer and/or per laboratory.

Offer valid for US and Canadian customers only. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. Offer valid while supplies last. Not valid on prior purchases. DWK distributors and employees are prohibited from participating. Void in government accounts or where prohibited by law or company policy. Prices subject to change without notice. Trademarks contained herein are the property of their respective owners.



DWK Life Sciences
1501 North 10th Street
Millville, NJ 08332
United States

US & Canada: 800.225.1437
Int'l: 856.825.1100
Fax: 856.825.1368
www.dwk.com