



**Buy 1 PIPET-PAL Controller & a Case of Serological Pipettes
 Receive a FREE case of DURAN® TILT Bottles**

Serological Pipets

Glass Serological Pipets

- Non-pyrogenic, plugged, sterile
- Manufactured from 51 expansion borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class B requirements.
- With negative graduations (Neg. Grad)

Plastic Serological Pipets

- Sterile, non-pyrogenic, and DNase/RNase-free. Irradiation sterilized for use with tissue culture work.
- Clean, plugged, sterile and ready for use.
- With negative graduations (Neg. Grad)

Catalog No.	Description	Case Qty
Sterile Glass		
72108-5110	Glass sero pipet, wide tip, 5 mL, TD	400
72108-10110	Glass sero pipet, wide tip, 10 mL, TD	400
72100-11100	Glass sero pipet, Yellow, 1 mL, Grad 0.01 mL, Neg. Grad 0.2 mL, TD, multi-pack	1000
72100-21100	Glass sero pipet, White, 2 mL, Grad 0.01 mL, Neg. Grad 0.2 mL, TD, multi-pack	500
72100-5110	Glass sero pipet, Blue, 5 mL, Grad 0.1 mL, Neg. Grad 1 mL, TD, multi-pack	500
72100-10110	Glass sero pipet, Orange, 10 mL, Grad 0.1 mL, Neg. Grad 2.5 mL, TD, multi-pack	500
Non-Sterile Glass		
72120-11100	Non-sterile, yellow, bulk pack, 1 mL, Grad 0.01 mL, Neg. Grad 0.2 mL, TD	1000
72120-21100	Non-sterile, white, bulk pack, 2 mL, Grad 0.01 mL, Neg. Grad 0.2 mL, TD	500
72120-5110	Non-sterile, blue, bulk pack, 5 mL, Grad 0.1 mL, Neg. Grad 1 mL, TD	500
72120-10110	Non-sterile, orange, bulk pack, 10 mL, Grad 0.1 mL, Neg. Grad 2.5 mL, TD	500
Plastic Sterile Pipets		
56900-11100	Plastic sero pipet, 1 mL, Grad 0.01 mL, Neg. Grad 0.3 mL, TD, individually wrapped	1000
56900-21100	Plastic sero pipet, 2 mL, Grad 0.01 mL, Neg. Grad 0.3 mL, TD, individually wrapped	500
56900-5110	Plastic sero pipet, 5 mL, Grad 0.1 mL, Neg. Grad 2 mL, TD, individually wrapped	200
56900-10110	Plastic sero pipet, 10 mL, Grad 0.1 mL, Neg. Grad 3 mL, TD, individually wrapped	200
56900-25210	Plastic sero pipet, 25 mL, Grad 0.2 mL, Neg. Grad 10 mL, TD, individually wrapped	200
56900-50510	Plastic sero pipet, 50 mL, Grad 0.5 mL, Neg. Grad 10 mL, TD, individually wrapped	100
56900-10010	Plastic sero pipet, 100 mL, Grad 1 mL, Neg. Grad 5 mL, TD, individually wrapped	50

PIPET-PAL™ Controller


- Lightweight (195g) and well-balanced, making even prolonged pipetting very comfortable
- Quiet, robust pump for high pipetting speed
- Precise finger-control of pipetting enables accurate dispensing, even drop-by-drop
- Versatile for both fine dosing and gravimetric dispensing
- Designed for use with all serological pipettes in the volume range 1-100mL
- Eliminate contamination with a replaceable, sterile filter
- Lithium-Polymer battery technology for long operating time
- Autoclavable pipette holder for cleanliness and safety
- Supplied with a universal power supply for worldwide use
- 1 year warranty

Catalog No.	Description	Case Qty
W155091	PIPET-PAL Controller, Translucent Blue	1

DURAN® TILT Media Bottle System



Cat. No.	Capacity (mL)	Thread GL	d (mm)	h (mm)	Case Qty
21 891 4459	500	56	55	151	4



Available at

Pipette.com

info@pipette.com • 800-243-3232

To claim your redemption (redeem by 7/31/19)

1. Obtain your Proof of Purchase (dated between 1/1/19 - 6/30/19). Packing slip, purchase order and invoice are all considered valid for proof of purchase.
2. Visit www.dwkdeals.com
3. Attach your Proof of Purchase and fill-out the online form.
4. Enter the part numbers for freeware. Free items must be of equal or lesser value. Visit www.DWK.com for more details. Not for resale. Limited to two redemptions per promotion, per customer and/or per laboratory.

Offer valid for US and Canadian customers only. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. Offer valid while supplies last. Not valid on prior purchases. DWK distributors and employees are prohibited from participating. Void in government accounts or where prohibited by law or company policy. Prices subject to change without notice. Trademarks contained herein are the property of their respective owners.



DWK Life Sciences
1501 North 10th Street
Millville, NJ 08332
United States

US & Canada: 800.225.1437
Int'l: 856.825.1100
Fax: 856.825.1368
www.dwk.com