

Most Popular

Pipette Service Order Form

Customer Information

FULL NAME

PHONE

E-MAIL

BILLING E-MAIL

SHIP TO ADDRESS

COMPANY NAME

DEPARTMENT

ADDRESS

CITY, STATE & ZIP CODE

BILL TO ADDRESS Check if same as "SHIP TO ADDRESS"

COMPANY NAME

DEPARTMENT

ADDRESS

CITY, STATE & ZIP CODE

PAYMENT INFORMATION

PO NUMBER OR CREDIT CARD NUMBER

EXP DATE / Master Card Visa Amex

***Please include a PDF copy of your PO**

QUOTE NUMBER

MY QUALER ACCOUNT (optional)

MY PIPETTE.COM ACCOUNT (optional)

Shipping/Handling Options

USE MY CARRIER AND MY ACCOUNT # FOR SHIPPING

Carrier Account #

FREE 2-DAY FEDEX SHIPPING WITH QUANTITIES OF: 15 PLUS FOR CALIBRATION OR 20 PLUS FOR VERIFICATION for smaller quantities, shipping rates below apply:

Ground Shipping (one way) \$15

2nd Day Shipping \$25

Overnight Shipping (one way) \$45
Continental USA (other areas call for quote)

Rates shown above are based on Domestic Shipping Rates Only

Select One Service Level	Verification*	Calibration	Calibration PLUS	Accredited
Manual Single Channel Pipette	\$27	\$39	\$44	\$54
Electronic Single, Repeater, Dispenser	\$28	\$44	\$49	\$57
Motorized Controller	\$28	\$44	N/A	N/A
6 Channel Pipette (manual/electronic)	\$46	\$46	\$51	\$118
8 Channel Pipette (manual/electronic)	\$46	\$67	\$72	\$149
12 Channel Pipette (manual/electronic)	\$46	\$78	\$83	\$181
16 Channel Pipette (manual/electronic)	\$46	\$88	\$93	\$212
24 Channel Pipette (manual/electronic)	\$46	\$99	\$104	\$244

Service Includes	Verification*	Calibration	Calibration PLUS	Accredited
Leakage, Functional & Calibration Check	Included	Included	Included	Included
Preventive Maintenance (PM) **	n/a	Included	Included	Included
As Found Readings	4x2	n/a	4x2	4x2
As Left Readings	n/a	4x2	4x2	10x2
Calibration Report***	Included	Included	Included	ISO Format
Comprehensive CAL Label	Included	included	Included	Included

Options	Verification*	Calibration	Calibration PLUS	Accredited
3rd Point Reading	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	<input type="checkbox"/> \$10	<input type="checkbox"/> \$12
Extended 10x2 Readings	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	<input type="checkbox"/> \$14	<input type="checkbox"/> \$7
Multi-Channel Readings (per channel)				
<i>As Found</i>			<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
<i>As Left</i>	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	Included
• <i>3rd Point Reading (per channel)</i>			<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
<i>As Found</i>			<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
<i>As Left</i>	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
• <i>Extended 10x2 Readings (per channel)</i>			<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
<i>As Found</i>			<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
<i>As Left</i>	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	Included
GMP/GLP Calibration Certificate	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	ISO Cert
DNA/RNA Decontamination (per pipette)	<input type="checkbox"/> \$15	<input type="checkbox"/> \$15	<input type="checkbox"/> \$15	<input type="checkbox"/> \$15

* If pipette fails Verification, maintenance service will be performed under Calibration Plus level.
 ** PM (Preventive Maintenance) includes new Seal & O-ring as specified by manufacturer, piston polishing, thorough inside & outside cleaning.
 *** Included with returned pipettes / Uncontrolled copy of GMP/GLP available online at no charge. The laboratory implements the shared risk approach, based on the binary simple acceptance decision rule, as defined in ILAC G8 when making compliance statements.

Calibration Cycle

MONTHS UNTIL NEXT CALIBRATION (Default is 6 Months)

3 Months 6 Months 12 Months Custom

PIPETTE QUANTITY

24 HOUR SERVICE (\$10 Per Pipette), limit of 25 pipettes per order, over 25 pipettes, please contact service@pipette.com
Transit time is not included, (business days only)

Notify me if repairs exceed \$60 per pipette. Pipettes that require rare parts may be delayed.

Paperless Account: Contact us for Qualer Access to Manage your pipettes and other assets for FREE!

*Calibration Warranty: 90 days labor / 1 year installed parts. Warranty not included with verification service.

REQUIRED

PIPETTES WERE DECONTAMINATED FROM RADIOACTIVE, CHEMICAL AND BIO HAZARDOUS MATERIALS WITH AN APPROPRIATE SOLUTION.

_____ CUSTOMER SIGNATURE


PLEASE PRINT THIS FORM AND INCLUDE EQUIPMENT LIST WITH YOUR SHIPMENT

PLEASE PRINT THIS FORM AND INCLUDE EQUIPMENT LIST WITH YOUR SHIPMENT

Pipettes For Trade- In:

Model Number	Serial Number	User/ID/Notes


Customer Notes and/or Requirements to Meet Your Internal SOP's are Available on Request



We have BIG news!

We have expanded and now have a calibration lab in Cleveland Ohio! This makes it easier and faster for us to service your pipettes and cuts down on shipping times to the eastern and central US!

Additional and customizable options such as multiple extra data points or readings are available. Please contact service@pipette.com or (800)243-3232 to expand your calibration service.



CONTACT US FOR QUALER ACCESS TO MANAGE YOUR PIPETTES AND OTHER ASSETS FOR FREE!