

Most Popular

Pipette Service Order Form

PLEASE PRINT THIS FORM AND INCLUDE WITH YOUR SHIPMENT

Select One Service Level	Verification*	Calibration	Calibration PLUS	Accredited Plus ISO 17025
Manual Single Channel Pipettes	\$26	\$37	\$42	\$51
Electronic Singles, Repeaters, Dispensers	\$27	\$42	\$47	\$54
Motorized Controller	\$27	\$42	N/A	N/A
6 Channel Pipettes (manual/electronic)	\$44	\$44	\$49	\$112
8 Channel Pipettes (manual/electronic)	\$44	\$64	\$69	\$142
12 Channel Pipettes (manual/electronic)	\$44	\$74	\$79	\$172
16 Channel Pipettes (manual/electronic)	\$44	\$84	\$89	\$202
24 Channel Pipettes (manual/electronic)	\$44	\$94	\$99	\$232

Service Includes	Verification*	Calibration	Calibration PLUS	Accredited Plus ISO 17025
Leakage, Functional & Calibration Check	Included	Included	Included	Included
Preventive Maintenance (PM) **	n/a	Included	Included	Included
As Found Readings	4x2	n/a	4x2	4x2
As Left Readings	n/a	4x2	4x2	10x2
Calibration Report***	Included	Included	Included	ISO Format
Comprehensive CAL Label	Included	Included	Included	Included

Options	Verification*	Calibration	Calibration PLUS	Accredited Plus ISO 17025	
3rd Point Reading	As Found	<input type="checkbox"/> \$5	n/a	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5
	As Left	n/a	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	<input type="checkbox"/> \$7
3rd Point Reading Per Channel	As Found	<input type="checkbox"/> \$5	n/a	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5
	As Left	n/a	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
Extended 10x2 Readings (per channel)	As Found	<input type="checkbox"/> \$7	n/a	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
	As Left	n/a	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	included
Multi-Channel Readings (per channel)	As Found	<input type="checkbox"/> \$7	n/a	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7
	As Left	n/a	<input type="checkbox"/> \$7	<input type="checkbox"/> \$7	included
DNA/RNA Decontamination (per pipette)	n/a	<input type="checkbox"/> \$15	<input type="checkbox"/> \$15	<input type="checkbox"/> \$15	
Printed copy of GMP/GLP Cal Cert.	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	<input type="checkbox"/> \$5	ISO Cert.	

*If pipette fails Verification, maintenance service will be performed under Calibration Plus Level.
 ** PM (Preventive Maintenance) includes new Seal & O-ring as specified by manufacturer, piston polishing, thorough inside & outside cleaning.
 *** Included with returned pipettes / Uncontrolled copy of GMP/GLP available online at no charge. The laboratory implements the shared risk approach, based on the binary simple acceptance decision rule, as defined in ILAC G8 when making compliance statements.

Calibration Cycle

MONTHS UNTIL NEXT CALIBRATION (Default is 6 Months)

3 Months 6 Months 12 Months Custom

PIPETTE QUANTITY

24 HOUR SERVICE (\$10 Per Pipette), limit of 25 pipettes per order, over 25 pipettes, please contact service@pipette.com
 Transit time is not included, (business days only).

Notify me if repair exceeds \$60 per pipette. Pipettes that require rare parts may be delayed.

Paperless Account: Contact Us for a Qualer Access to Manage your Pipettes and other Assets for FREE!
 *Calibration Warranty: 90 days labor / 1 year installed parts.

PIPETTES WERE DECONTAMINATED FROM RADIOACTIVE, CHEMICAL AND BIO HAZARDOUS MATERIALS WITH AN APPROPRIATE SOLUTION

_____ CUSTOMER SIGNATURE

LIST TRADE-IN PIPETTES AND SPECIAL REQUIREMENTS ON REVERSE SIDE
PLEASE PRINT THIS FORM AND INCLUDE WITH YOUR SHIPMENT

Customer Information

FULL NAME _____

PHONE _____

E-MAIL _____

BILLING E-MAIL _____

SHIP TO ADDRESS

COMPANY NAME _____

DEPARTMENT _____

ADDRESS _____

CITY, STATE & ZIP CODE _____

BILL TO ADDRESS Check if same as "SHIP TO ADDRESS"

COMPANY NAME _____

DEPARTMENT _____

ADDRESS _____

CITY, STATE & ZIP CODE _____

PAYMENT INFORMATION

PO* NUMBER OR CREDIT CARD NUMBER

_____|_____|_____|_____|_____|_____|_____|_____|_____|_____|_____|_____|_____|_____|_____|_____|

EXP DATE

____|____|/____|____| Master Card Visa Amex

MY QUALER ACCOUNT (optional) _____

MY PIPETTE.COM ACCOUNT (optional) _____

PROMO CODE (optional) _____

*Please include a PDF copy of your PO

Shipping/Handling Options

USE MY CARRIER AND MY ACCOUNT # FOR SHIPPING

Carrier _____ Account # _____

FREE 2-DAY FEDEX SHIPPING WITH (Domestic shipping only):
 5 PLUS PIPETTES FOR CALIBRATION OR 10 PLUS FOR VERIFICATION
 For smaller quantities, domestic shipping is \$11 each way

Overnight Shipping (each way) _____ \$28
 Continental USA (other areas call for quote)

PLEASE PRINT THIS FORM AND INCLUDE WITH YOUR SHIPMENT

List Pipettes Below or Contact Us for Qualer Access to Manage Your Pipettes and Other Assets for FREE!

Model Number	Serial Number	User/ID/Notes

Pipettes for Trade-In: Do you have any pipettes for trade-in? Please fill out the table below.

Customer Notes and/or Requirements to Meet Your Internal SOP's are Available on Request

Additional and customizable options such as multiple extra data points or readings are available.
Please contact service@pipette.com or (800)243-3232 to expand your calibration service.

