

#### Detachment

Remove the 4 safety pins by slightly rotating and pulling by the closed end.

Remove the draw latches by pulling up-out & off of the dewar shell latch catches.

The dewar shell and the Sliderite™ pallet can then be detached if required.

Reverse process to reattach.



**24-HOUR HELP LINE**  
**(949) 470-2305**

Cryoport, Inc. • 17305 Daimler St. • Irvine, CA 92614  
Phone: 1-949-470-2300 Fax: 1-949-470-2306  
www.cryoport.com • E-mail Contact: cs@cryoport.com