

WHITNEY

Try
Valitome™!

The newest way to work with bone cement

whitnecurettes.com
800.338.4237

The Valitome's design was inspired by orthopedic surgeon,
Dr. Wayne Goldstein, Illinois Bone and Joint Institute.
To show our appreciation we named our new product
after his granddaughter, Valerie.

WHITNEY VALITOME™

WHITNEY
Has A New Tip For You!



The newest way to work with bone cement

Osteotome-style Beveled Tip

Straight edge applies cement evenly, thoroughly for better adherence.
Narrow profile gets into hard-to-reach places.
Sharp edge removes cement easily, completely, safely from flat surfaces.

Made of High Strength Plastic

Minimizes worry about scratching surface of prosthesis.

Ergonomic Shape

Ease of handling facilitates precise movements.
Molded-in thumb indentation instantly orients tip of tool.

Color Coded

Immediate identification of proper tool eliminates confusion.

Disposable

Saves valuable time – no need to clean cement from tool.

Arrives Sterile

Individually packaged, ready to go.

NO-SCRATCH™ WHITNEY CURETTES™ – Available in Four Styles



Whitney Curette™ - Small
Order # OR1201
Tip size: 8 x 6 mm
Color: Blue
Quantity: Carton of 25



Whitney Sculps Knife™
Order # OR1203
Tip Size: 20 x 14 mm
Color: Green
Quantity: Carton of 25



Whitney Curette™ - Large
Order # OR1202
Tip Size: 20 x 12 mm
Color: Yellow
Quantity: Carton of 25



Whitney Valitome™
Order # OR1205
Tip Size: 13mm Osteotome
Color: Teal
Quantity: Carton of 2

*The distinctive look of Whitney Surgical Tools constitutes protected trade dress under U.S. and International law.

Valitome™



Real Solutions for Real Problems™