

# BIOHAZARD KEEPERS™ – FLOOR MODEL

Manage Your Biohazards  
With WHITNEY MEDICAL SOLUTIONS:  
Real Solutions for Real Problems™



You can collect and dispose of biohazardous lab waste conveniently and economically in the floor model Biohazard Keeper™ with the unique wide-mouth Reclose-A/Remove-A-Lid™ cover.

## COLLECT AND DISPOSE OF LAB WASTE EASILY AND SAFELY

The Terminal® brand Biohazard Keeper™ offers an easy and economical way to collect and dispose of lab waste. The sturdy corrugated paperboard construction offers strength, and the 2 mil red plastic liner provides leak-resistance. Dispose of non-sharp biohazardous lab waste, including petri dishes, tubes, plates, pipets and pipet tips, specimen cups, and tubing sets.

The unique Reclose-A/Remove-A-Lid™ allows you to reclose the container during use to reduce odors and mess, or you can remove the lid entirely. The wide-mouth opening allows you to dispose of large items, or entire racks of tubes or pipet tips at once. When the container is full, a separate flap can be folded over the top to permanently close the container for final disposal.

Containers store flat, conserving your storage space, and assemble easily in seconds. Warnings are printed on the carton in English and Spanish.

CAUTION: Not a sharps container—do not use for disposal of contaminated sharps.

*biohazards*  
You'll wonder how you ever managed <sup>1</sup> without us.™

## TERMINAL® brand BIOHAZARD KEEPER™

**Dimensions:** 12" W x 12" D x 27" H

**Capacity:** 63 liters

**Shipping Weight:** 14 lbs.

**Quantity:** Case of 6

**Order Number:** BH2009

**Features:** Sturdy corrugated construction, heavy-duty 2 mil liner, and unique wide-mouth opening with Reclose-A/Remove-A-Lid™ cover.



Real Solutions for Real Problems™