

BIOHAZARD BENCH-TOP KEEPERS™

Manage Your Biohazards
With WHITNEY MEDICAL SOLUTIONS:
Real Solutions for Real Problems™



You can get rid of biohazard waste right at your workstation without making unnecessary trips to your central biohazard waste container with the plastic-lined Bench-Top Keeper™.

COLLECT WORKSTATION WASTE ECONOMICALLY

Easily collect contaminated waste items like petri dishes, pipets, test tubes, applicator sticks, swabs and other non-sharp items right at your workstation using Terminal® brand Biohazard Bench-Top Keepers™. Less expensive than rigid sharps containers, Bench-Top Keepers save you money and keep your workstation neat and clean.

Bench-Top Keepers have a plastic lining and hold a moderate amount of liquid. Self-locking tops keep contents safely enclosed. They store flat – requiring far less space than rigid plastic containers – then fold into shape for use, with interlocking flaps and orange biohazard symbols. Autoclave contaminated materials in the container, then dispose of according to your local, state or federal regulations. Available in two sizes.

You'll wonder how you ever managed ^{biohazards} without us.™

LARGE (shown)

Dimensions: 6"L x 6"W x 10"H

Capacity: 3000cc

Shipping Weight: 15 lbs.

Quantity: Case of 50

Order Number: BH2003

Features: Free Standing; webbed corners and plastic lining keep liquids from soaking through.

SMALL

Dimensions: 4"L x 4"W x 7"H

Capacity: 1400cc

Shipping Weight: 12.5 lbs.

Quantity: Case of 100

Order Number: BH2002

Features: Same great features as large size; takes up less bench space.



Real Solutions for Real Problems™