

BIOHAZARD PIPET KEEPER™

Manage Your Biohazards
With WHITNEY MEDICAL SOLUTIONS:
Real Solutions for Real Problems™



You can dispose of pipets and other long objects up to 15 inches in length easily with the Pipet Keeper™.

DISPOSE OF LONG ITEMS EASILY

Use our biohazard Pipet Keeper™ to collect contaminated pipets, swabs, and other long objects up to 15" in length right at your workstation. Terminal Brand Pipet Keepers are made of sturdy paper-board printed with orange biohazard symbols. They are lined with plastic to avoid soak-through and designed with webbed corners to hold moderate amounts of liquid. Priced lower than rigid plastic sharps containers, they are disposable, saving time and making it easier to maintain high lab standards. An interlocking top with reclosable flaps makes the Pipet Keeper easy to use whether you are right or left handed. The Pipet Keepers are shipped flat for efficient storage and can easily be folded together for use.

biohazards
You'll wonder how you ever managed ¹without us.™

Terminal® brand Pipet Keeper™

Dimensions: 3" W x 3" H x 15½" L

Capacity: 2250 cc

Shipping weight: 8 lbs.

Quantity: Case of 50

Order number: BH2006

Features: Collects contaminated long items right at your workstation.



Real Solutions for Real Problems™