

SAFE-KEEPER™

Manage Your Biohazards
With WHITNEY MEDICAL SOLUTIONS:
Real Solutions for Real Problems™



You can save valuable space on your bench top with the hangable and leak-resistant Safe-Keeper™

SAVE SPACE AND MONEY WITH OUR SAFE-KEEPER™

Use our Safe-Keeper™ containers to collect biohazardous waste right at your workstation without cluttering up your bench top. You can dispose of capillary tubes, pipet tips, swabs, stir rods, and other contaminated items safely and economically in a Terminal® brand Safe-Keeper. Made of sturdy plastic-lined paperboard, the Safe-Keeper has a unique solid bottom design, unlike other look-alike products. The solid bottom keeps contents from leaking or falling out. Safe-Keepers can be autoclaved or incinerated when filled.

Safe-Keepers are shipped flat to save space and set up easily for use. Each case comes with a Hook-Up™ self adhesive hooks so you can hang your Safe-Keepers on the wall or a shelf, keeping your bench top clear. Additional hooks sold separately; suitable for keys, lab coats, goggles, etc.

You'll wonder how you ever managed ^{biohazards} without us.™

SAFE-KEEPER™

- Dimensions:** 5"L x 2"W x 11"H
- Capacity:** 1100cc
- Shipping Weight:** 13 lbs.
- Quantity:** Case of 200
- Order Number:** BH2005

HOOK-UP™

- Dimensions:** 2"L x 2"W x 3/4"H
- Shipping Weight:** 1 lbs.
- Quantity:** Box of 4
- Order Number:** BH6012
- Features:** Great for keys, light items, hanging extra Safe-Keepers.



Features: Solid bottom construction avoids leaks; can be wall-mounted.



Real Solutions for Real Problems™