

## Measuring/Indications

PRODUCT	INDICATIONS	SIZING INFORMATION
1	<ul> <li>OA</li> <li>Meniscal tears</li> <li>Tendinopathy</li> <li>Residual pain/inflammation (post-op)</li> </ul>	Circumference of knee at patella w/ knee at 90 degrees
Y	<ul> <li>Post-surgical acute swelling for TJR and arthroscopy:         <ul> <li>5-7 days post, patient wear 4-6 weeks</li> </ul> </li> <li>Post-op swelling that hasn't responded to PT and/or hosiery, 4-6 weeks post</li> </ul>	Circumference of patient thigh at widest point
	<ul><li>Plantar fasciitis</li><li>Tendinopathy</li></ul>	Based on patient's shoe size; please indicate gender, i.e. "8/F"
	<ul> <li>Lateral and medial epicondylitis</li> <li>Residual pain/inflammation post UCL reconstruction</li> </ul>	Circumference at elbow
	<ul> <li>Distal swelling from shoulder trauma</li> <li>Post-surgical acute inflammation</li> </ul>	Circumference at mid-bicep
	Carpal tunnel	Circumference at wrist
I	<ul> <li>Undersleeve for any walking boot application - coded as L2840*</li> <li>*not payable in all jurisdictions</li> </ul>	Shoe size - reference sizing guide re: calf circumference
1	<ul> <li>Suspension sleeve for any functional knee brace - coded as L2397</li> </ul>	Circumference of thigh - 4" above patella center

Fax: 530-267-4882