



Yes, you can.®

1900 Bathing Seat Lift



1900S Shown

The Invacare 1900 Bathing Seat Lift is specifically designed for use with 3652G and 3602GXL Side-Entry Bathing Tubs. Safely transfer residents with dignity and minimal handling. The soft seat is warm to the touch and offers true comfort.


1900 Bathing Seat Lift



Convenience
An optional digital scale is accurate to nearest pound or kilogram.



Maneuverability
Hospital-grade casters provide easy maneuverability, lockable casters are standard equipment.



Ease of use
Vertical motion is operated by rugged, ergonomic hand pendant.



Peace of mind
Each unit comes with **two** 24V DC sealed gel-cell batteries and convenient wall mounted charger.

Specifications

- Lifting capacity: 350 lbs.
- Seat height: 18"
- Vertical range of motion: 23"
- Batteries: (2) 24V sealed gel-cell
- Base length: 47"
- Base width: 28.5"
- Maximum height: 73"

Features

- Resident support arms rotate up to 120° for easy transfers
- Heavy gauge steel mobile base is powder coated for durability and aesthetics
- Contoured, soft, closed-cell polyurethane foam seat with perineal opening designed for resident comfort and easy cleaning
- Ergonomically designed pushbar
- Powerful and quiet DC lift actuator drive mechanism

Invacare Continuing Care, Inc.

www.invacare-cc.com

info_icc@invacare.com

USA

One Invacare Way
Elyria, Ohio 44035
USA
800-668-2337

Canada

994 Hargrieve Road
London, Ontario
N6E 1P5 Canada
800-668-2337

Trademarks are identified by the symbol™ and ®.
All trademarks are owned or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

© 2010 Invacare Corporation. Form 04.12.5/10



Yes, you can: