


TheraPure™ 6302G & 6302GADH Supine Tubs



Connection	Location	Type
C	10" (254 mm) from centerline, 5.5" (140 mm) up from floor	3/4" Cold water stub (with recommended shut off valve)
D	5" (127 mm) from centerline, 29" (736 mm) from tub back edge	2" Drain (floor drain recommended)
E	10.5" (267 mm) from centerline	115V AC, 60 Hz 15 A Input
H	6.5" (165 mm) from centerline, 5.5" (140 mm) up from floor	3/4" Hot water stub (with recommended shut off valve)

Invacare Continuing Care, Inc.

www.invacare-cc.com

info_icc@invacare.com

USA

One Invacare Way
Elyria, Ohio 44035
USA
800-668-2337

Canada

994 Hargrieve Road
London, Ontario
N6E 1P5 Canada
800-668-2337

Trademarks are identified by the symbol™ and®.
All trademarks are owned or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

© 2010 Invacare Corporation.



Yes, you can.