

3point products Design Line™ Thumb Arthritis Splint



The **DESIGN LINE THUMB ARTHRITIS SPLINT** combines comfort and fashion in a whole new way. Designed to support the CMC joint for relief of the grinding pain of CMC joint arthritis, it's easy-to-apply and fully adjustable for a customized fit. Available in 3 fashion patterns, this moderate support splint provides **COMFORT, FIT** and **FASHION**.

Features

- Contoured strap supports and carefully abducts the CMC joint to reduce the pain from CMC joint arthritis or ligament injury.
- Soft foam-lined stretch material contours to the hand with natural compression so there is no need to overtighten the splint.
- Straps can be secured anywhere on the hook-engageable material to achieve a custom fit.
- Three sizes and three patterns offer an ideal range of fit and fashion and are the perfect alternative to bulky, black-on-black splints.
- 100% breathable, 100% latex-free; hand or machine washable.

Order Information

Measure across the palm at the level of the big knuckles (MP joints).

Pattern	Size	Left	Right	Fits
Zebra	Small	P2008-L2ZB	P2008-R2ZB	2 3/8" to 3" (6cm-7.6cm)
Zebra	Medium	P2008-L3ZB	P2008-R3ZB	3" to 3 5/8" (7.6cm-9.2cm)
Zebra	Large	P2008-L4ZB	P2008-R4ZB	3 5/8" to 4 1/2" (9.2cm-11.4cm)
Leopard	Small	P2008-L2LP	P2008-R2LP	2 3/8" to 3" (6cm-7.6cm)
Leopard	Medium	P2008-L3LP	P2008-R3LP	3" to 3 5/8" (7.6cm-9.2cm)
Leopard	Large	P2008-L4LP	P2008-R4LP	3 5/8" to 4 1/2" (9.2cm-11.4cm)
Red Trim	Small	P2008-L2RB	P2008-R2RB	2 3/8" to 3" (6cm-7.6cm)
Red Trim	Medium	P2008-L3RB	P2008-R3RB	3" to 3 5/8" (7.6cm-9.2cm)
Red Trim	Large	P2008-L4RB	P2008-R4RB	3 5/8" to 4 1/2" (9.2cm-11.4cm)

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com