



3pp® Wrist P.O.P.™ Splint

The 3pp® Wrist P.O.P.™ (Point of Pressure) Splint applies adjustable compression over the distal ulna and a counterpoint of pressure under the distal radius for stability and effective treatment of TFCC injuries, ulnar-sided wrist pain, DRUJ (distal radioulnar joint) instability and arthritis.




Features

- Realigns and stabilizes the DRUJ to reduce painful symptoms of TFCC injuries, DRUJ.
- The corrective strap and pressure pads allow the user to position force across the distal ulna and radius to relieve painful clicking and carpal instability.
- Two 1/8" thick and two 1/4" thick pressure pads allow the user to choose the combination of pads needed for pain relief.
- Foam-lined base sleeve will not slip or roll and eliminates the need for sticky closures touching the skin.
- The Design Line cover adds a fashionable touch while preventing hook closures from catching clothing.
- Machine washable. Latex free.

Order Information

Measure around your wrist to size. Splint fits right or left arm.



	Product Number	Fits
Sm/Med	P3030-23	5" to 9"
Med/Lrg	P3030-34	9" to 11"

Fits right or left arm

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com