

# 3point: products®

Designed for Comfort, Fit and Function



**We've made the BEST even BETTER.**


Oval-8 splints are now easier on your eyes and easier on your patient's finger.



Look for the **plus +** sign on one end – that end makes the splint fit looser when put on on the finger first.



Improved molding makes the **size easier to read** on the bottom of the band.



**Smoother, seamless bands** feel even better on patient's fingers.



## Oval-8® Finger Splints

Oval-8 splints treat mallet fingers, swan neck, trigger finger, hypermobility and more. The open 3-point pressure design allow the splints to treat multiple conditions and can be custom fit in seconds.

### PACKAGING

Available in Single and Mixed size packages and in clinic Kits for fast, convenient fitting. Individual sizes available in packages of 1 or 5.

### Oval-8® Sizing Set



One of each of the 14 sizes. Use the Sizing Set to quickly and accurately find the size that fits securely and comfortably on the finger.

### Oval-8® Kit



The Oval-8 Kit includes a total of 44 splints and keeps the splints organized and easy to see what sizes you need to restock.


### Oval-8® Graduated Sets



Packaged with 3 consecutive sizes, Graduated Sets help assure your patient has the correct size, even if swelling or temperature affects size.

## 3pp® Step Splints™

Static progressive 3pp Step Splints treat IP deviation, IP flexion and IP extension contractures. Easy-to-fit and easy for your patient to adjust, make these a great go-to choice for corrective splinting.



3pp® Side Step Splint

Adjustable tension corrects IP joint deviation caused by arthritic deformity or injury. Easily adjusted with the non-slip strap to restore ROM.



3pp® Step Down Splint

Restore composite DIP and PIP flexion following injury or surgery. Allows wearer to apply stretch to tolerance for optimal compliance.



3pp® Step Up Splint

Restore PIP and DIP extension to reduce contractures caused by boutonniere deformities or tendon injuries.

All 3pp Step Splints are available in 3 sizes. Not recommended for the thumb.

## 3pp® Soft Finger Splints



3pp® Buddy Loops®

Simple loop and wrap design make these the preferred strap to protect & restore ROM.

Available in 3 widths:  
½", ¾" and 1"



3pp® Final Flexion Wrap™

Restore end range composite IP flexion after injury or dynamic splinting. Easy adjustment assures compliance & success.

One size fits most



3pp® Finger Trapper™

Provide comfortable, non-slip traction for dynamic splinting. Conforms around the nail without discomfort.


One size fits most

## CMCcareThumb Brace

**Basko** Healthcare

The CMCcareThumb Brace contours around the first metacarpal to stabilize the CMC joint for pain-free motion. Ideal for treating Type II swan neck deformities of the thumb, the CMCcare restores proper CMC alignment for improved function. Carefully designed to avoid pressure in the web space, and smooth rounded edges are softly flared for comfort.

Available in 3 sizes in right and left models.



Carefully molded to stabilize the CMC joint, the CMCcare is a superior solution to treat CMC arthritis.




Squeeze the embedded stay to contour it to stabilize the base of the 1st metacarpal.



The CMCcare comes with a replacement strap and interchangeable closure tabs. Addition straps are available.



Strap is removable and replaceable. Splint is packaged with 2 straps.



Use Gel Mate® silicone gel to cushion or relieve prominent joints.



Edges can be easily adjusted using a spot heat gun.

### Superior Design and Materials Supreme Comfort

#### Fix Comfort Wrist Brace Manutec® by Orliman

Luxuriously cushioned with pressure absorbent foam and covered with smooth, breathable velour cover and straps for the ultimate in comfort, Malleable aluminum is used to create a carefully structured brace that conforms to the contours of the hand and forearm.

Fully adjustable strapping system and closure tabs easily secure the splint in place for firm immobilization of the wrist with no restriction of finger movement.


Ambidextrous design fits both right and left hand



Add the malleable thumb attachment (sold separately) to convert to a Wrist and Thumb splint



Palmar contours allow full finger and thumb motion



Soft edges won't impinge in the web and the thumb is malleable for a custom fit

#### Fix Comfort Thumb Brace Manutec® by Orliman

Cushioned, breathable, velour covered splint. Malleable interior frame allows for a fully custom contoured brace. Provides firm MP and CMC support with the utmost comfort.

Fully cushioned lightweight splint provides firms CMC and MP control without limiting finger or thumb IP motion. Comfortably protects the thumb during daily tasks.

Ambidextrous design fits both right and left hand and comes in two sizes.





## 3pp® Thumb Splints

3pp Thumb Splints offer varying levels of control to reduce thumb pain from CMC arthritis, ligament injury, post-op protection and more.



3pp® ThumSling®

Contoured strap applies light compression at the CMC joint to relieve pain. Hook receptive material allows straps to attach anywhere for optimal fit.

Two sizes in a choice of black or gray



3pp® ThumSling® NP

The same great ThumSling design with a greater degree of support. Stretch material maintains support as the hand moves. Straps close anywhere on the splint for optimal fit.

Two sizes in black



3pp® ThumSpica Plus

The best combination of a custom molded orthosis and easy-to-fit, wrap-on support. Precut thermoplastic creates a custom-molded splint in seconds.

One size fits most

## 3pp™ Design Line™ Thumb Splint



Wrap-around strap lifts and supports the CMC joint to relieve the grinding pain of basal joint arthritis. Binding tape prevents edges from raveling; adjustable thumb strap eliminates compression at the IP; and a removable stay provides adjustable support at the MCP. The contoured fit of the Design Line Thumb Splint and memory stretch material provide support without bulkiness.

Three sizes- Small Medium & Large, make it easy to stock and they're available in 4 fashion options.

4 fashion options



Camo



Red Trim



Leopard



Zebra

### 3pp® Wrist P.O.P.™ Splint and 3pp® Carpal Lift

Choose the best approach for treating TFCC injury or Ulnar instability

#### 3pp® Wrist P.O.P. Splint



Use to apply adjustable counterforce pressure on the distal radius and ulna to stabilize the DRUJ. Includes camo cover to keep strap in place.

2 sizes fit right or left

#### 3pp® Carpal Lift



Use to boost the ulnar carpus to reduce instability without limiting motion. Palmar cushion relieves ulnar-sided wrist pain and Guyon's Canal Syndrome.

2 sizes in right and left models

### 3pp Ez FIT ThumSpica Splint

Treats:

De Quervain's, Tendinitis, Thumb Sprains

The 3pp Ez FIT ThumSpica Splint stabilizes the CMC and MP joints. Fits right or left hand.

- Easy one-handed application
- Softly cushioned with no hard edges
- One size fits most
- Adjustable stay is easily contoured



Fully Cushioned



Contoured fit



Malleable stay

## 3pp® Elbow P.O.P.™ Splint



The non-slip, easy-to-apply 3pp Elbow P.O.P. Splint applies an adjustable point of pressure to reduce pain from medial or lateral tendinosis without the need to over-tighten the splint to keep it in place. Relief without constriction.

Two sizes fit right or left.



Includes a black camo cover to keep straps in place and add a fashionable touch

## Ulnar Deviation Splints

### Polycentric Hinged Ulnar Deviation Splint



Designed to align the MP joints to reduce ulnar deviation that can be corrected with mild to moderate force. The dorsal based splint leaves the palm free and is hinged to allow full flexion and extension for performance of daily tasks.

4 sizes in right and left models

### Radial Hinged Ulnar Deviation Splint



The Radial Hinged Ulnar Deviation Splint adjusts each finger individually to improve the alignment of strongly deviated fingers. The radial hinge allows free metacarpophalangeal (MP) flexion and extension for improved hand function.

3 sizes in right and left models

### Comforter Splint



Softly padded, the Comforter Splint keeps the wrist supported and the MP joints aligned in neutral to reduce inflammation and morning stiffness caused by rheumatoid arthritis. The IP's and thumb are free to manage bedding and nighttime activities.

3 sizes in right and left models



### A Great Group of Products for Scar Care

**Gel Mate® Silicone Gel** medical grade silicone gel sheets are self-adherent on one side with an anti-microbial fabric backing that prevents clothing from sticking and pulling it away from the scar.

The **SkinSational® Massage Brush** safely and gently massages to reduce or prevent adhesions, improve mobility and aid in healing. The soft, pliable head is safe to use on recently healed scars.

The **Scar Rx Kit** combines Gel Mate, the SkinSational and the highest quality organic oil and lotion for a complete kit for successful scar care.



GelMate



Scar Rx™



SkinSational™

### Topricin® Pain Relief Products

No Odor No Sting No Irritation Just Pain Relief

Safe, effective non-prescription pain relief cream. Safe on all skin types, simply rub in the cream and experience relief from muscle strain, arthritis, nerve pain, foot pain and fibromyalgia.

Topricin works with no stinging, no odor, no stickiness, no smelly hands after application... And IT WORKS.

Specialized formulas for muscle & joint, foot or Fibromyalgia pain.





## 3pp® Toe Loops®

The non-slip loop and wrap design makes 3pp Toe Loops the ideal choice to treat hammertoes, fractures or post-op applications. Soft foam lining cushions for comfort. Wide Toe Loops fits the great toe. Narrow Toe Loops fit smaller toes.

Packages of 5 and 25



## 3pp® Bunion-Aider™

The 3pp Bunion-Aider corrects alignment of the great toe by applying a gentle stretch to stabilize the toe while relieving Hallux Valgus pain. Includes a removable stay for stronger correction.

One size fits most



## 3pp® Arch Lift™

The 3pp Arch Lift is designed to lift and support the arch to relieve pain from plantar fasciitis, flat foot, heel spurs and metatarsalgia.

One size fits most



## FS6 Foot Compression Sleeve

Comfortable, form-fitting knitted foot sleeve employs gradient compression to support the foot and ankle. Provides relief for pain from plantar fasciitis, Achilles tendinitis and heel spurs.

4 sizes - 2 sleeves per pack

## Heathcare Professionals



3-Point Products offers healthcare customers professional wholesale pricing on all our products. Our professional pricing program include:

- Professional wholesale pricing up to 30% off advertised retail pricing
- Quantity discounts on many products
- Free shipping on orders above \$500
- New product sampling discounts

Contact customer service to set up your account today to take advantage of SAVINGS

**Call 888-378-7763**

**Mon-Fri 8:30 am to 5:00 pm EST**

**[service@3pointproducts.com](mailto:service@3pointproducts.com)**

For complete ordering & product information, visit  
[3pointproducts.com/Resources](http://3pointproducts.com/Resources)

---

## Patient Direct Purchasing



We work hard to earn your trust and provide you and your patients exemplary service and quality products. Our Oh My Arthritis division was created as a consumer website your patients can trust, and where they can easily find the products you have recommended.

We have gone one step further and created HCP referral coupons, so you can help your patients save on the products you recommend.

### PATIENT SAVING COUPONS

Simply provide your patients with a coupon from our coupon pads or let them know to use the code HCP15 for their first order.



Facsimile coupon - not for patient use

*Thank you*  
for viewing our product catalog

888-378-7763

Mon-Fri 8:30 am – 5:00 pm (EST)

[service@3pointproducts.com](mailto:service@3pointproducts.com)

For details on our full line of products, visit

[www.3pointproducts.com](http://www.3pointproducts.com)