



Gel Mate Silicone Gel Sheet



Gel Mate is a medical grade silicone gel sheet for scar management. Gel Mate is self-adherent on one side and backed by an anti-microbial (bacteria resistant) cloth that glides under clothing on the other. Helps flatten and soften raised scars and may prevent hypertrophic (thickened) or Keloid scars. Silicone gel sheets are 100% latex-free

Features

- Medical grade silicone gel helps improve the color and appearance of new or mature scars.
- Applies long-term light pressure on scars and has been show to add a hydrating effect to the scar and surrounding skin that helps in increasing mobility and faster healing.
- Reduce itching and discomfort experienced when scars begin to heal.
- Sheets include a measurement grid on the back to determine the correct size.

Order Information

Size	Product Number
4" x 6"	P9510 (1 sheet)
4" x 6"	P9510-5 (5 sheets)
2.5" x 5"	P9511 (1 sheet)



Watch our instructional video on the Scar Rx Kit to learn more about Gel Mate:
www.3pointproducts.com/product-videos-of-splints/scar-rx-video



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com