



## Gel Tubes



Vitamin enriched Gel Tubes ease pressure and friction on fingers and toes. Use Gel Tubes to protect and cushion corns, hammertoes, finger lacerations, ulcerations, enlarged joints and more.

### Features

- Protect, cushion and help secure splints in place when worn under custom-made or prefabricated splints like Oval-8 Finger Splints.
- Superior durability, these washable and reusable tubes allow for repeated wear.
- Hypoallergenic oil based polymer releases medical grade mineral oil and vitamins to soften and protect skin and assist with scar remodeling.
- Tubes are 1/2" (1.5 cm) wide, 1 5/8" (4 cm) long and stretch well.
- Hand washable and reusable.
- Product Number P1111-5 (pack of 5). One size fits most.



**Call 888-378-7763**

Mon-Fri 8:30am - 5:00pm EST for ordering details

**[www.3pointproducts.com](http://www.3pointproducts.com)**