



Mueller Green Fitted Wrist Brace



The Mueller Green Fitted Wrist Brace comfortably cushions and supports weak or injured wrists while maintaining full range of motion of the thumb and fingers. Padded dorsal and volar stays combine with two adjustable straps to provide firm support to hold the wrist in a neutral position to help relieve the pain and swelling associated with Carpal Tunnel Syndrome or ligament injuries.

Features


- Supports and stabilizes the wrist with 2 adjustable straps and stays for firm support.
- This easy-to-fit and easy-to wear orthosis is comfortable for all day and/or all night wear for relief from Carpal Tunnel Syndrome, ligament injuries and wrist sprains.
- Breathable fabric is protected with AEGIS Microbe Shield Controls to protect against odors, stains, and product wear and tear.
- Environmentally friendly, this orthosis uses recycled and non-petroleum materials, and is latex-free.
- Two sizes for the right and for the left hand minimize clinic stocking.
- Firm Support

Order Information

Measure wrist circumference to size

	Right	Left	Fits
Small/Medium	P3232-R23	P3232-L23	5"-8" (12.7-20cm)
Medium/Large	P3232-R45	P3232-L45	8"-10" (20.3-25.4cm)

Watch our **instructional video on the Mueller Green Fitted Wrist Brace** on our website:
<http://www.3pointproducts.com/product-videos-of-splints/mueller-green-fitted-wrist-brace-video>



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com