



Mueller Night Support Wrist Brace



The Mueller Night Support Wrist Brace holds the wrist in a neutral position during sleep. A beaded cushion in the palm limits fisting and the circumferential design provides dorsal and volar support for optimal relief of Carpal Tunnel Syndrome, arthritis and tendinitis.

Features

- 2 metal stays and 360 degree support provide rigid stability to hold the wrist in neutral.
- Beaded cushion supports the palm to limit fisting for optimal relief of tingling and pain from Carpal Tunnel Syndrome.
- Nighttime support to relieve swelling and stiffness caused by arthritis or tendinitis.
- Fully cushioned for optimal comfort during sleep
- Antimicrobial materials help keep the splint clean. Latex-free materials
- One Size Fits Most – can be worn on the left or right hand
- Firm Support

Order Information

Measure wrist circumference to size

One size fits wrist circumference from 5 ¾" to 9" (14.6 cm to 22 cm)

Universal design fits right or left hand

Product Number: P3233-U

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com